

0963139-000000

## of

Kirk M. Ririe,

and

for

## CONTAINER FOR CARRYING OUT AND MONITORING BIOLOGICAL PROCESSES

Attorney Docket 7475-66667

Attorneys:

Steven R. Lammert  
Jill T. Powlick  
BARNES & THORNBURG  
11 South Meridian Street  
Indianapolis, Indiana 46204

# CONTAINER FOR CARRYING OUT AND MONITORING BIOLOGICAL PROCESSES

## BACKGROUND

5

### 1. Related Applications.

A This is a divisional of U.S. Patent Application Serial No. 08/869,275,  
filed June 4, 1997, entitled System And Method For Monitoring PCR Processes which  
is a continuation-in-part of U.S. Patent Application Serial No. 08/658,993, filed June  
10 4, 1996 (now abandoned) entitled System And Method For Monitoring PCR  
Processes which is a continuation-in-part of U.S. Patent Application Serial No.  
08/537,612, filed October 2, 1995, entitled Method For Rapid Thermal Cycling of  
Biological Samples which is a continuation-in-part of U.S. Patent Application Serial  
No. 08/179,969, filed January 10, 1994 (now U.S. Patent No. 5,455,175), entitled  
15 Rapid Thermal Cycling Device which is a continuation-in-part of U.S. Patent  
Application Serial No. 07/815,966 filed January 2, 1992 (now abandoned) entitled  
Rapid Thermal Cycling Device which is a continuation-in-part of U.S. Patent  
Application Serial No. 07/534,029 filed June 4, 1990 (now abandoned), entitled  
Automated Polymerase Chain Reaction Device, each of the above-identified  
20 applications are now each individually incorporated herein by reference in their  
entireties.

The copending U.S. application filed in the U.S. Patent and Trademark  
on June 4, 1997 entitled Monitoring Hybridization During PCR as serial no.  
08/869,276 and naming Carl T. Wittwer, Kirk M. Ririe, and Randy P. Rasmussen as  
25 inventors is hereby incorporated by reference in its entirety.

### 2. The Field of the Invention.

This invention relates generally to apparatus which are used to carry  
out biological processes, such as the polymerase chain reaction. More specifically,  
30 the present invention relates to apparatus and methods which carry out thermal cycling  
and monitoring of various biological reactions, such as the polymerase chain reaction.

### 3. The Background Art.

In numerous areas of industry, technology, and research there is a need to reliably and reproducibly subject samples to thermal cycling. The need to subject a sample to repeated thermal cycles is particularly acute in biotechnology applications.

- 5 In the biotechnology field, it is often desirable to repeatedly heat and cool small samples of materials over a short period of time. One such biological process that is regularly carried out is cyclic DNA amplification.

Cyclic DNA amplification, using a thermostable DNA polymerase, allows automated amplification of primer specific DNA, widely known as the  
10 "polymerase chain reaction" or "PCR." Automation of this process requires controlled and precise thermal cycling of reaction mixtures usually contained in a plurality of containers. In the past, the container of preference has been a standard, plastic microfuge tube.

- Commercial programmable metal heat blocks have been used in the  
15 past to effect the temperature cycling of samples in microfuge tubes through the desired temperature versus time profile. However, the inability to quickly and accurately adjust the temperature of the heat blocks through a large temperature range over a short time period, has rendered the use of heat block type devices undesirable as a heat control system when carrying out processes such as the polymerase chain  
20 reaction.

- Moreover, the microfuge tubes which are generally used have disadvantages. The material of the microfuge tubes, their wall thickness, and the geometry of microfuge tubes is a hindrance to rapid heating and cooling of the sample contained therein. The plastic material and the thickness of the wall of microfuge  
25 tubes act as an insulator between the sample contained therein and the surrounding medium thus hindering transfer of thermal energy. Also, the geometry of the microfuge tube presents a small surface area to whatever medium is being used to transfer thermal energy. The continued use of microfuge tubes in the art, with their suboptimal geometry, indicates that the benefits of improved thermal transfer (which  
30 come by increasing the surface area of a sample container for a sample of constant volume) has heretofore not been recognized.

Furthermore, devices using water baths with fluidic switching, (or mechanical transfer) have also been used as a thermal cyclers for the polymerase chain reaction. Although water baths have been used in cycling a polymerase chain reaction mixture through a desired temperature versus time profile necessary for the reaction to take place, the high thermal mass of the water (and the low thermal conductivity of plastic microfuge tubes), has been significantly limiting as far as performance of the apparatus and the specificity of the reaction are concerned.

Devices using water baths are limited in their performance. This is because the water's thermal mass significantly restricts the maximum temperature versus time gradient which can be achieved thereby. Also, the water bath apparatus has been found to be very cumbersome due to the size and number of water carrying hoses and external temperature controlling devices for the water. Further the need for excessive periodic maintenance and inspection of the water fittings for the purpose of detecting leaks in a water bath apparatus is tedious and time consuming. Finally, it is difficult with the water bath apparatus to control the temperature in the sample tubes with the desired accuracy.

U.S. Patent No. 3,616,264 to Ray shows a thermal forced air apparatus for cycling air to heat or cool biological samples to a constant temperature. Although the Ray device is somewhat effective in maintaining a constant temperature within an air chamber, it does not address the need for rapidly adjusting the temperature in a cyclical manner according to a temperature versus time profile such as is required for biological procedures such as the polymerase chain reaction.

U.S. Patent No. 4,420,679 to Howe and U.S. Patent No. 4,286,456 to Sisti et al. both disclose gas chromatographic ovens. The devices disclosed in the Howe and Sisti et al. patents are suited for carrying out gas chromatography procedures but do not provide thermal cycling which is substantially any more rapid than that provided by any of the earlier described devices. Rapid thermal cycling is useful for carrying out many procedures. Devices such as those described in the Howe and Sisti et al. patents are not suitable for efficiently and rapidly carrying out such reactions.

In particular, the polymerase chain reaction (PCR) is a fundamental DNA amplification technique essential to modern molecular biology. Despite its

00631329 000300



usefulness and popularity, the current understanding of PCR is not highly advanced. Amplifications must be optimized by trial and error and protocols are often followed blindly. The limited understanding of PCR found in the art is a good example of how those skilled in the art are content to utilize a powerful technique without reflection or  
5 comprehension.

Biological processes such as PCR require temperature cycling of the sample. Not only does the prior art, as explained above, carry out temperature cycling slowly, the prior art also ignores the underlying principles which allow PCR to work and could be used to make PCR even more useful. Thus, it would be a great advance  
10 in the art to provide methods and apparatus which are particularly adaptable for rapidly carrying out PCR and analyzing the reaction which is taking place, particularly if such reaction is analyzed as it is taking place, that is, in real time.

#### BRIEF SUMMARY AND OBJECTS OF THE INVENTION

15 In view of the above described state of the art, the present invention seeks to realize the following objects and advantages.

It is an object of the present invention to provide an apparatus for accurately controlling the temperature of biological samples.

It is a further object of the present invention to provide a thermal  
20 cycling apparatus for quickly and accurately varying the temperature of biological samples according to a predetermined temperature versus time profile.

It is another object of the present invention to provide an apparatus suitable for subjecting a number of different biological samples to rapid thermal cycling.

25 It is also an object of the present invention to provide a thermal cycling apparatus having a thermal transfer medium of low thermal mass which can effectively subject samples to a large temperature gradient over a very short period of time.

It is a further object of the present invention to provide an apparatus  
30 which can subject a biological sample to rapid thermal cycling using air as a thermal transfer medium.

5

10

15

20

25

30

09639-0600

In order to rapidly cool the sample chamber, the preferred apparatus includes a means for forcing air into the sample chamber and a means for dispersing the air forced into the sample chamber. The preferred structures included in some embodiments are a high velocity fan which functions to force air into the sample chamber and a rotating paddle which functions to disperse the air into the chamber. In some embodiments, a means for venting allows the air to escape from the sample chamber taking the unwanted heat with it. The present invention allows heating and cooling of a sample to take place both quickly and uniformly.

In accordance with the method and the apparatus of the present invention, a control structure provides means for operating the system through a desired time versus temperature profile. The present invention is particularly well suited for carrying out automated polymerase chain reaction procedures.

The controller of the present invention allows the biological samples to pass through a predetermined temperature cycle corresponding to the denaturation, annealing and elongation steps in the polymerase chain reaction. In use, the apparatus of the present invention allows rapid optimization of denaturation, annealing, and elongation steps in terms of time and temperature, and shortened time periods (ramp times) between the temperatures at each step.

The present invention particularly decreases the total time required for completion of polymerase chain reaction cycling over prior art thermal cycling devices while at the same time significantly increasing specificity and yield.

In accordance with another aspect of the present invention, the present invention provides methods and apparatus for monitoring of DNA amplification so as to track the progress of such procedures. In particular, the present invention provides methods and apparatus for continuous fluorescence monitoring of the polymerase chain reaction procedure. In preferred embodiments of the present invention, optical components are combined with structures to provide rapid temperature cycling in order to continuously monitor DNA amplification by a variety of different fluorescence techniques. Glass capillary sample containers and composite plastic/glass sample containers allow rapid heat transfer from the preferred thermal transfer medium (allowing 30 amplification cycles in less than 15 minutes when a gas

such as air is used as the thermal transfer medium) and simultaneous monitoring of the reaction.

In accordance with another aspect of the present invention, optical techniques are used to monitor the progress of the reaction as the reaction is ongoing.

5 In some preferred embodiments of the invention, fluorescent probes are added to the reaction mixture. The present invention then monitors the fluorescence at least once during a temperature transition, and preferably the fluorescence is acquired two or more times during a temperature transition, either from a single sample or from multiple samples. In some preferred embodiments a rotating carousel is included to  
10 sequentially move the samples, one-by-one, to a monitoring location with all of the samples being simultaneously subjected to rapid thermal cycling. Desirably, embodiments of the present invention provide for monitoring of fluorescence once per amplification cycle or monitoring temperature, time, and fluorescence continuously throughout each amplification cycle.

15 Using the present invention, a 3-dimensional plot of temperature, time, and fluorescence, can be obtained. Fluorescence vs. temperature plots of hybridization probes discriminate between the cumulative, irreversible signal of exonuclease cleavage and the temperature-dependent, reversible hybridization of adjacent probes. Hybridization probes are more useful than hydrolysis probes because  
20 the temperature dependence of fluorescence can be followed and used to detect alterations in product sequence, *i.e.*, polymorphisms and mutations. Using dyes that fluoresce in the presence of double stranded DNA, product denaturation, reannealing and extension can be followed within each cycle. The present invention provides apparatus and methods for rapidly carrying out DNA amplification reactions which  
25 combines amplification and analysis of the reaction in under fifteen minutes and more preferably in under fifteen minutes and most preferably in under ten minutes.

#### BRIEF DESCRIPTION OF THE DRAWINGS

In order to better appreciate how the above-recited and other  
30 advantages and objects of the invention are obtained, a more particular description of the invention briefly described above will be rendered by reference to specific embodiments thereof which are illustrated in the appended drawings. Understanding

that these drawings depict only typical embodiments of the invention and are not therefore to be considered limiting of its scope, the invention will be described and explained with additional specificity and detail through the use of the accompanying drawings in which:

Figure 1 shows a perspective view of a thermal cycling apparatus adapted for thermal cycling of biological samples and adapted especially for use in cyclic DNA amplification, according to the concepts of the present invention.

Figure 2 is a side elevation view of the fluid chamber portion of the apparatus of Figure 1.

Figure 3 is an interior plan view of the fluid chamber portion of the apparatus illustrated in Figure 1.

Figure 4 shows an interior plan view of the fluid chamber of another embodiment of the present invention.

Figure 5 shows an optimized temperature versus time profile for a  
15 polymerase chain reaction using the thermal cycling device of the present invention.

Figure 6 shows graphically the effect of denaturation time on polymerase chain reaction yields using one thermal cycling device of the present invention.

Figure 7 shows graphically the effect of annealing time on polymerase  
20 chain reaction specificity and yields using the thermal cycling device of the present  
invention.

Figures 8A-B, which are perspective and elevational cross sectioned views, respectively, of another preferred embodiment of the present invention.

Figure 8C is a diagrammatic representation of the relationship of the  
25 heat producing element and the capillary tubes holding the biological samples in the  
embodiment illustrated in Figures 8A-B.

Figure 9A shows the results of four different temperature/time profiles (A-D) and their resultant amplification products after thirty cycles (A-D).

Figure 9B shows a cycle of another preferred temperature/time profile  
30 used by the present invention.

Figures 9C-G show exemplary cycles of other preferred temperature/time profiles used by the present invention.

Figure 10 provides a block diagram of a temperature slope control circuit in accordance with the present invention.

Figure 10A is a graphical representation of the effect of the temperature transition rate from the product denaturation temperature to the primer  
5 annealing temperature on reaction product specificity.

Figure 11 is a schematic view of a preferred rapid temperature cycler with fluorescence detection in accordance with the present invention.

Figure 11A is a temperature v. time chart of showing one preferred operation of the apparatus of Figure 11.

10 Figure 12 is a representation of three dimensional plots of temperature, time, and fluorescence during amplification of a hepatitis B DNA fragment in the presence of SYBR Green I.

Figures 12A-C are representations of two dimensional plots of temperature vs. time, fluorescence vs. time, and fluorescence vs. temperature which  
15 are together shown as a three dimensional plot in Figure 12.

Figure 13 is a fluorescence vs. temperature projection during the amplification of a 536 base pair fragment of the human  $\beta$ -globin gene in the presence of SYBR Green I.

Figure 14 is a cycle number vs. fluorescence plot obtained in  
20 accordance with an aspect of the present invention.

Figure 14A provides a legend for Figure 14, and subsequent figures, indicating different initial template copy numbers.

Figure 15 is a cycle number vs. fluorescence ratio plot obtained in accordance with an aspect of the present invention.

25 Figure 16 is a fluorescence ratio vs. temperature plot obtained in accordance with one aspect of the present invention.

Figure 17 is a fluorescence ratio vs. temperature plot obtained in accordance with one aspect of the present invention.

Figure 18A is a graph representing an equilibrium PCR paradigm.

30 Figure 18B is a graph representing a kinetic PCR paradigm.

Figure 18C is a graph representing different time/temperature profiles near an annealing temperature.

0063139-080300

Figures 19A-19D are representations of different sample container configurations.

5

preferred sample container in accordance with the present invention.

10

Figure 20 shows the optical layout of another preferred embodiment in accordance with the present invention to provide continuous monitoring of a sample undergoing DNA amplification.

15

Figures 21A-D show composite plastic/glass containers into which biological samples are loaded.

20

Figure 22A charts the effectiveness of light piping by viewing the tip rather than the side of capillary sample container.

Figure 22C is a high level block diagram showing the tasks which are performed by one preferred embodiment of the present invention which includes a rapid temperature cycler with epifluorescence detection.

Figure 22E is a plot of fluorescence vs. time for a PCR reaction in which fluorescence feedback was used to control reaction parameters.

Figure 23 is a plot of fluorescence vs. time showing showing the inverse relationship between temperature and fluorescence.

Figure 24 is a plot of temperature vs. time showing the inverse relationship between temperature and fluorescence.

5                Figure 25 is a plot of fluorescence vs. temperature for three different PCR products in the presence of SYBR Green 1 acquired during a 0.2 degree per second temperature transition through the product melting temperatures.

Figure 26 is a plot of fluorescence vs. time showing product annealing for different concentrations of PCR product in the presence of SYBR Green 1.

10              Figures 27A and 27B are cross sectional schematic views of the embodiment represented in Figure 28 in a run mode and a load mode ,respectively.

Figure 28 is a schematic representation of another embodiment of the present invention which is a rapid temperature cycler with fluorescence detection at the tip of the sample containers and which includes positioning for fluorescence  
15              detection in two dimensions to optimize detection.

Figure 29 is a perspective view of the exterior of the embodiment of the present invention including the components illustrated in the schematic representation of Figure 28.

20              Figures 30A-30V are detailed schematic diagrams of the electrical components of one preferred embodiment of the present invention.

Figures 31A and 31B are perspective and cross sectional views, respectively, of a sample handling system in accordance with the present invention.

Figure 32 is a schematic representation of another embodiment of the present invention which accommodates multiple sample handling trays.

25

#### DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENTS

Reference will now be made to the drawings wherein like structures will be provided with like reference designations.

30              As shown in Figure 1, the one preferred thermal cycling device 10 includes a closed loop fluid (most preferably air) chamber, generally designated at 11, which is adapted to accept samples to be cycled through vent door 14. The closed loop fluid chamber 11 includes a plurality of compartments each of which will be



described shortly. The device 10 also includes a controller 12 which can be programmed by means of input keys 25 and display 26 to cause the chamber 11 to be cycled through a series of temperatures over a predetermined period of time. The thermal cycling of chamber 11 can be used to carry out numerous procedures and is particularly suited for amplification of primer specific DNA from samples containing reaction mixtures as will be explained below.

The closed loop fluid chamber 11 is enclosed in a generally box shaped configuration by housing 13. Blower mounting boards 16, if desired, can be located so as to section off a smaller rectangular section of the chamber 11 and function to support and secure a generally cylindrically shaped lower housing 15 thereto. Alternatively, the fan of the blower 28 may be housed integrally within chamber housing 13.

The interior of blower housing 15 contains the blades and shaft of the blower. The blower motor (not shown) is located externally of blower housing 15, and therefore exteriorly of the enclosed chamber 11. In this configuration, the blades and shaft are the only parts of the blower which become exposed to the circulating hot fluid within chamber 11. It would be disadvantageous to mount the motor within the chamber which would subject the motor to temperature variations and also would add the thermal mass of the motor to that which is subject to heating and cooling. The reduction of thermal mass exposed to the fluid in chamber 11 is desirable to the overall performance of the device 10 in its function of subjecting samples placed therein to a desired temperature versus time profiles, using either predetermined profiles or by altering one or more reaction parameters as the reaction continues, as will be more fully explained below.

The blower 28 is a well known type of blower usually identified as an "in line" type blower which preferably employs a propeller type fan, due to its generally low thermal mass, or if desired, a squirrel cage type fan, the fan preferably having a 75 cubic feet per minute minimum capacity.

The solenoid platform 17 has secured thereto a solenoid 18. The solenoid armature 19 is attached to upper end 21 of rod 20 which is rigidly attached to vent door 14 and rotatably attached to housing 13 at points above and below the vent

door 14. The rod 20 therefore allows vent door 14 to freely rotate relative to the housing 13 about the rod's longitudinal axis.

A spring 22 is attached at one of its ends to the housing 13 by support post 23. The opposite end of spring 22 is attached to the top end 21 of rod 20 directly adjacent the attachment of solenoid armature 19. The spring 22 is drawn between these two attachment points so as to be in tension. The spring 22 therefore tends to draw top end 21 toward the support post 23, which in turn tends to rotate vent door 14 to its closed position. When solenoid 18 is actuated, armature 19 tends to pull top end 21 of the rod 20 in the direction of the solenoid 18, which is opposite the direction of pull of spring 22, and which tends to open the vent door 14.

Controller, generally designated at 12, is electrically attached to the chamber 11 by means of a transmission cable 24. The cable 24 also supplies power to the blower motor (not shown), and to the heat coil 31. Further, the controller 12 also is connected to thermocouple sensor 35 for receiving signals corresponding to temperature data, and to solenoid 18 for triggering the solenoid armature.

Controller 12 can be any well known type of temperature controller unit which is programmable to control the heat coil 31, vent door 14, and blower so as to achieve predetermined temperatures as a function of time within the chamber 11, and which is also capable of being programmed to actuate a relay output for driving a solenoid at predetermined time periods and chamber temperature levels. A preferred temperature controller 12 for use in the embodiment of Figures 1-3 is a Partlow MIC-6000 proportional temperature controller, available through Omega Engineering Inc, of Stanford, Connecticut, as the Model No. CN8600 process controller.

As shown in Figures 2 and 3, the interior of chamber 11 is sectioned off into four main compartments. The blower compartment 28 is formed of the blower housing 15 and the blower mounting plates 16. The entirety of blower compartment 28 is filled with the fan and shaft portions of a blower as has been described above. The blower can be any of a number of well-known designs, as has been described above, and has therefore been omitted from Figure 3 for purposes of clarity. It is sufficient for the present invention to understand that the fan located in blower compartment 28 draws fluid into the blower compartment 28 through inlet opening 36 and pushes the fluid out of exit opening 37.

It is preferred that the fluid be driven by the blower at a rate of at least 75 cubic feet per minute. It is important however, in regard to the present invention, to realize that the fluid located in chamber 11 only contacts the fan and a portion of the drive shaft of the blower, the blower motor itself being located outside of the blower housing 15 so as to avoid any contact thereof with fluid in the chamber 11. This consideration contributes to the speed of operation of the invention to minimize the material which contacts the fluid inside the chamber 11 so as to minimize the thermal mass of material which must be heated and/or cooled thereby during the cycling process. By minimizing the thermal mass which must be heated or cooled by the fluid, the response time necessary to bring the contents of chamber 11 to a uniform temperature is greatly diminished.

Fluid exiting blower compartment 28 through outlet opening 37 enters heating compartment 29. Fluid passing into heating compartment 29 must pass by heating coils 31. If the heating coils 31 get hotter than the fluid passing into heating compartment 29, the fluid will become heated thereby as it is forced through the compartment. The heating coil is preferably a 1,000 watt (125 VAC) nichrome wire coil wound around a microsupport. However, any heating unit suitable for heating the type of fluid present in the chamber may be used. The particular heating coil of embodiment of Figures 1-3 is manufactured by Johnstone Supply, of Portland, Oregon.

The heating coil is activated by an output relay included in the controller 12. The preferred relay is a 25 A, 125 VAC solid state relay manufactured by Omega Engineering Inc. of Stamford, Connecticut as Model No. Omega SSR 240 D25.

Fluid passing through heating compartment 29 becomes incident on baffles 32 and 33 before passing into the reaction compartment 30. Baffles 32 and 33 tend to break up any laminar fluid flow and generate turbulence therein to effectively mix the fluid so that it arrives in reaction compartment 30 at an homogenous temperature.

Thermocouple sensor 35 provides an electrical input signal to controller 12 which corresponds to the fluid temperature in the reaction compartment 30. Temperature monitoring during operation of the thermal cycling device 10 is

5                   The fluid passing from the reaction compartment 30 to the return air  
compartment 34 must pass through sample compartment 27 (as shown in dashed  
lines). Sample compartment 27 will also be explained momentarily.

The fluid in return compartment 34 has been slightly cooled due to the heat transfer therefrom into samples located in sample compartment 27. The fluid in return compartment 34 is drawn through inlet opening 36 into blower compartment 28 where it is again forced, by action of the fan, out through outlet opening 37 into the heating compartment 39. Thus, the fluid chamber 11, when operating with vent door 14 closed, is a closed loop fluid chamber which continuously recirculates the fluid along a closed loop path through each compartment thereof in order to bring the contents therein to a uniform temperature. Continuous circulation of the air in the air chamber 11 allows the samples in sample compartment 27 to be brought to a predetermined temperature as quickly as possible, and then to be held at that temperature, if desired.

When the device 10 must be used to not only heat material located in the reaction compartment 27, but also to subsequently cool these materials as quickly as possible to a temperature at or above the ambient fluid (air) temperature, the controller 12 can be programmed to actuate solenoid 18 to cause vent door 14 to open and allow large quantities of ambient fluid to immediately flood the compartment 11 while heated fluid therein simultaneously escapes.

25 Deactivation of the heating coil 31 while continuing activation of the  
blower with vent door 14 open, will draw ambient fluid into return compartment 34  
and from there into the blower compartment 28. The blower will then push this  
ambient fluid through heating compartment 29 where it will pass directly into reaction  
compartment 30 without being heated by coil 31. The ambient fluid then passes  
30 through the sample compartment 27 and escapes out of chamber 11 through the vent  
door 14. Due to the minimum thermal mass of material located in chamber 11, and  
the action of the blower fan, vast quantities of ambient fluid will be forced past the

sample compartment 27, and from there out of the chamber 11. Thus, rapid cooling of samples or material located in the reaction compartment 27 is obtained.

The sample compartment 27 is sized so as to allow a plurality of samples, such as hollow elongate glass tubes containing a sample therein, to be easily located in a spaced apart orientation so that fluid may be evenly distributed around each sample. If desired, the sample compartment 27 may be sized and configured so as to allow insertion of a rack, basket, or the like which has been configured so as to accept a plurality of samples in uniform spaced apart configuration so as to simplify loading the samples into the sample chamber 27.

Access to sample compartment 27 is accomplished by rotation of the vent door 14 to its open position. Once the vent door 14 is rotated to approximately 90 degrees from its closed position, the sample compartment 27 is easily accessible there through. Also, as can be seen in Figures 1-3, rotation of vent door 14 approximately 90 degrees from its closed position causes return fluid compartment 34 to be substantially closed off from the reaction compartment 30. Thus, when the device 10 of the present invention is in a "cooling" mode, ambient fluid enters directly into the return fluid compartment 34 and is forced through the blower compartment 28, heating compartment 29, reaction compartment 30, and sample compartment 27 substantially along the same path as the closed loop fluid flow path described above. The fluid is then forced out of the air chamber 11 and prevented from passing back into air return compartment 34 by the positioning of the vent door 14 between the sample compartment 27 and the return fluid compartment 34.

Thus, the vent door 14 not only allows ambient fluid to enter the chamber 11, it can also prevent the fluid from recirculating in a loop fashion through the chamber 11. Instead, fluid is forced to pass through the sample compartment 27 and then out of the chamber 11 to aid in the rapid cooling of the sample contents and chamber 11.

When the device 10 of the present invention is used for cyclic DNA amplification, repetitive cycling through different temperatures is required. Samples containing a reaction mixture for the polymerase chain reaction generally must be cycled approximately 30 times through temperature changes which correspond to the denaturation, annealing and elongation phases of the amplification process.

The device 10 of the present invention, due to its novel characteristics described above, is capable of cycling samples in significantly shortened periods compared to the prior art. For example, the DNA amplification application of the embodiment represented in the figures can pass through a temperature versus time  
5 profile cycle in 30-60 seconds (see Figure 5). This same cycle using prior art devices would take approximately 5-10 times longer. These low cycle times have proven also to increase yield and specificity of the polymerase chain reaction over prior art cycling.

10

---

### Example 1

The polymerase chain reaction was run in a 10  $\mu$ l volume with 50 ng of human genomic template DNAs, 0.5 mM of each deoxynucleotide, 500 nM of each of two oligonucleotide primers GGTTGGCCAATCTACTCCCAGG (SEQ ID NO: 5) and GCTCACTCAGTGTGGCAAAG (SEQ ID NO: 6) in a reaction buffer  
15 consisting of 50 mM Tris-HCl (pH 8.5 at 25°C), 3.0 mM magnesium chloride, 20 mM KCl, and 500  $\mu$ g/ml bovine serum albumin. Thermus aquatics DNA polymerase (0.4  $\mu$ ) was added, the samples placed in 8 cm long, thin-walled capillary tubes (manufactured by Kimble, Kimax 46485-1), and the ends fused with a laboratory gas burner so that an air bubble was present on both ends of each tube.

20

The capillary tubes were then placed vertically in a holder constructed of 1 mm thick "prepunched perfboard" (manufactured by Radio Shack). The mixture was cycled 30 times through denaturation (90-92°C), annealing (50-55°C), and elongation (72-75°C) to give the temperature versus time profile of Figure 5. Temperature monitoring of the capillary tubes was done with a miniature  
25 thermocouple (IT-23, Sensortek, Clifton, NJ) placed in 10  $\mu$ l of deionized water and connected to a thermocouple monitor (BAT-12, Sensortek). Amplification products were fractionated by electrophoresis on a 1.5% agarose gel. Specific amplification products were obtained in good yield.

---

30

Due to the fact that the device 10 of the present invention uses air as the thermal transfer medium instead of water, it has the advantage that heat transfer occurs through a low heat capacity medium (air) which can be warmed very rapidly.

The response time for sample cooling is very fast due to the use of thin walled glass capillary tubes for holding samples, instead of plastic microfuge tubes as has been done in the past with prior art processes, and by minimizing the thermal mass of material inside the chamber 11 (see Figure 5). Such response times can allow for optimization of the time and temperature requirements for the denaturation, annealing, and elongation steps in the polymerase chain reaction.

Further, shortened "ramp" times are obtained, i.e., the time required to bring the temperature of the sample from one temperature level to the next temperature level corresponding to phases in the amplification process is shortened. This decreases the time required for a complete amplification, as well as allowing specific study of annealing, denaturation and enzyme kinetics within a polymerase chain reaction protocol.

The baffles 32 and 33 (as shown in Figure 3) may be used if desired to achieve better temperature homogeneity within the sample compartment 27. As shown in this embodiment, baffles 32 and 33 decrease the temperature variation in the reaction compartment 30 from about 10°C, to about 2°C. If desired, further (or more complicated) baffles may be used to further decrease the temperature variation in reaction compartment 30. Alternately, as shown in Figure 4 the fan may be positioned downstream from the heating coil 31, but before the sample compartment 27 to achieve more uniform mixing.

Amplification products obtained through the use of apparatus 10 are at least qualitatively and quantitatively as desirable as those obtained through the manual water bath cycling method. However, advantages in specificity and yield are possible with rapid thermal control of the reaction mixture.

Figure 6 shows the effect of the temperature versus time profile of Figure 5 as used with the thermal cycling apparatus 10 on specificity (i.e., one specific product yield as opposed to a plurality of similar or "shadow" products). As can be seen, the shorter the ramp and annealing time, the greater the product specificity. The

rapid temperature response of the apparatus 10 allows improved specificity and yield which is not possible with prior art systems.

Figure 7 shows the effect of varying the denaturation time of the temperature versus time profile of Figure 5 as used with the thermal cycling apparatus 10 of the present invention on DNA amplification yields. The brighter vertical lines each correspond to a particular time at a denaturation temperature. As can be seen, the yield is greatest at the shortest possible denaturation time. Such a result is not possible with prior art systems.

As has been shown, by decreasing the thermal capacity (thermal mass) of the apparatus 10, the present invention can markedly decrease the total time required for carrying out the polymerase chain reaction. In addition, the use of small sample volumes further shortens the total time required for the reaction and also reduces the amounts of expensive reagents which must be used by up to about 90%, thus further reducing the cost of carrying out procedures using the present invention. For example, in the embodiment represented in Figures 1-3, capillary tubes 108 having inner diameters in the range from about 0.25mm to about 1.0mm can desirably be used. In some applications, capillary tubes 108 having inner diameters in the range from about 0.02mm to about 0.1mm can also be desirably used.

The apparatus 10 of the present invention is useful for amplifying DNA from any source. Although particular configurations and arrangements of the present invention have been discussed in connection with the specific embodiments of the thermal cycling device 10 as constructed in accordance with the teachings of the present invention, other arrangements and configurations may be utilized. For example, various fluids other than air, of generally low thermal mass, may alternatively be used in the device 10.

Another embodiment of the present invention is represented in Figures 8A-C. Figure 8A is a perspective view and Figure 8B is an elevational cross sectioned view of the additional embodiment. It will be understood that many of the earlier explained components and teachings also have application in the embodiment illustrated in Figures 8A-C. Thus, only the pertinent additional information concerning this embodiment will be provided below. Importantly, in the embodiment of Figures 8A-C, the heat producing element is adjacent to the biological sample



containers allowing faster heating and cooling of biological samples as explained below.

As will be appreciated shortly, the apparatus of Figures 8A-C provides even greater improvement over the prior art in the speed at which thermal cycling can be carried out, e.g., 15 or 30 cycles of DNA amplification in 30, 15, 10, or even fewer, minutes. Furthermore, the apparatus 100 provides better thermal homogenization throughout the samples than previously possible.

Shown in Figure 8A is the general configuration of the housing 102 of the embodiment. The housing 102 rests on feet 104 (best seen in Figure 8B) and functions to hold the other described structures in place and to isolate those structures which become hot from the surrounding environment. Included in the embodiment 100 of Figure 8A are input keys 25 and a display 26 as in the previously described apparatus 10. The previously described control structures can readily be modified or used as a pattern for a control means for use in the embodiment of Figures 8A-C.

As shown best in the cross sectional view of Figure 8B, a sample chamber is designated by bracket 106. A lid 138 connected to the housing 102 by a hinge 131 can be opened to allow access to the sample chamber 106. The sample chamber 106 is preferably cylindrical in shape but can be of any shape or size required by the particular application.

The sample chamber 106 is preferably lined with a black colored foam material 110 whose surface has light absorbing characteristics with the bulk of the thickness of the foam having insulating characteristics. The black foam material can be one which is readily available in the art and one fabricated from a plastic material. The foam 110 is preferably a material which is readily cooled by the air passing there over, i.e., the material has low thermal conductivity and a porous surface.

The dark or black porous surface of the material converts shorter wavelength radiation striking the surface into longer wavelength radiation, i.e., heat, which is radiated into the sample chamber.

The foam 110 functions to thermally isolate the sample chamber from the surrounding air space in the housing and also to convert the light emitted by lamp 112 into thermal energy. The foam 110 can be replaced with other structures. For example, a material having a black, dark, or other nonreflective surface, such as a thin

sheet of polycarbonate having one surface painted black, can be backed by an insulative material, such as a fiberglass or foam material. The black or dark surface, which can be painted on a number of different substrates, converts shorter wavelength radiation striking it into thermal radiation while the insulative material thermally isolates the sample chamber from the surrounding environment. Thus, using the teachings provided herein, those skilled in the art can utilize many different materials and structures as a lining for the sample chamber.

The lamp 112 is preferably a 500 watt halogen lamp. If appropriate control devices are used, higher power lamps or a plurality of lamps, such as four 500 watt halogen lamps, can be used. A lamp socket 112A is attached to the housing 102 by a support 112B. The lamp 112 is able to very rapidly and uniformly heat the sample chamber 106 to the desired temperature. Other sources of heat, i.e. infrared radiation, such as the earlier described nichrome wire element, can also be used within the scope of the present invention.

Represented in Figure 8B are two thin-walled capillary tubes 108 such as those described earlier. While two thin-walled capillary tubes 108 are shown, the sample chamber 106 can hold many such tubes. The thin-walled capillary tubes 108 have several important advantages over previously used devices as described earlier and, together with the sample chamber 106, function as the one presently preferred example of a means for holding a biological sample.

It will be appreciated that many other structures performing equivalent or similar functions can also be used. The thin-walled capillary tubes 108 are preferably left partially extending out of the sample chamber through apertures 140 for ease of access but may be completely contained within the sample chamber 106 as may numerous other fluid holding structures which are suited to particular applications. The preferred thin-walled capillary tubes 108 have a capacity of about  $10\ \mu\ell$ . As will be understood, the volume of the sample should be kept small, and the surface area of the sample holding structure relatively large, and together they present a relatively small thermal mass. It is also preferred that the sample holding structure contain a volume anywhere from about  $1\ \mu\ell$  to about  $10,000\ \mu\ell$  but those skilled in the art will appreciate that other volumes of samples can also be used within the scope of the present invention if the different thermal mass of the structure is considered.

The lamp 112 and the insulative foam 110 together provide rapid and uniform heating of the sample contained in the thin-walled capillary tubes 108 and the air contained within the sample chamber 106. A thermocouple 134 is included within the sample chamber 106 to sense the temperature within the chamber and is used to maintain the desired temperature within the sample chamber as earlier described.

The thermocouple 134 is preferably one available in the art whose thermal response substantially matches the thermal response of the biological sample and the container holding the same. Such thermocouples can be commercially obtained from sources such as Idaho Labs which manufactures a line of thermocouples referred to as metal sheathed, J-type thermocouples. The matching of the thermal response of the thermocouple to that of the biological sample and container can be preferably carried out by inserting a micro thermocouple, such as the model IT-23 thermocouple available from PhysiTemp as known in the art, into a typical biological sample being held by the chosen container and subjecting the sample and the thermocouple under test to the same temperature changes. The thermocouple under test, or some external criteria, can be changed until the thermal response of the thermocouple suitably matches the thermal response of the sample and its container.

The arrangement represented in Figure 8B provides more uniform heating and cooling of the sample than previously available devices. In previously available devices, transfer of heat throughout the sample is carried out by convection through the sample. Convection induced movement of the sample within whatever structure is used to hold the sample is caused by temperature gradients or differences in the generally small biological samples (e.g., 10-100  $\mu\ell$ ).

The effect of temperature gradients within the sample become more pronounced and more difficult to control as the cycle time for a sample decreases. The existence of uneven temperatures within a sample, and particularly the reliance on "mixing by convection" within the sample relied upon by the prior art devices, generally increases the cycle time for a sample and likely has deleterious effects on the biological sample. The apparatus 100 is capable of providing heating and cooling such that thermal differences within a 10  $\mu\ell$  sample are maintained at not greater than  $\pm 1^\circ\text{C}$  at all times during a 30 second cycle.

In order to promote uniform heating and cooling, it is preferred that the thin-walled capillary tubes 108 be at least somewhat uniformly spaced from the heat source, for example, lamp 112 in apparatus 100. Figure 8C provides a diagrammatic top view of the lamp 112 and the plurality of thin-walled capillary tubes 108 as arranged in the apparatus 100 represented in Figures 8A-B.

In the arrangement represented in Figure 8C, the thin-walled capillary tubes 108 which are farthest from the lamp 112 (as indicated by line F) are preferably no more than substantially 40%, and more preferably no more than substantially 25%, farther from the lamp 112 than the distance between the lamp 112 and those thin-walled capillary tubes 108 which are closest to the lamp 112 (as indicated by line N). For example, the distance indicated by line N can be about 7.3 cm while the distance indicated by line F can be about 8.5 cm.

It will be appreciated that the arrangement of the thin-walled capillary tubes 108 can be other than that represented in the figures, for example, circular or semi-circular. Moreover, it will be appreciated that the point from which to measure the distance between the heat producing element and the sample containers will vary as the type and size of the heat producing element varies. For example, the heat producing element may comprise a plurality of lamps or electric resistive elements which vary in shape and size. In some embodiments, it may also become important to consider the distance from the sample chamber wall the sample containers are positioned. In the illustrated embodiment, the apertures 140 (see Figure 8A) function as a means for holding the sample containers but other structures performing equivalent functions can also be used in accordance with the present invention.

The apparatus 100 also cools the samples contained in the capillary tubes 108 very rapidly and uniformly. In order to cool the sample chamber 106, air from outside the housing 102 is drawn into the interior of the housing through a lower housing portal 114 by a fan 116 which is connected to a motor shaft 122 driven by a motor 118. Since rapid cooling of the sample chamber is desired, it is preferred that the combination of the motor 118 and the fan 116 be able to move sufficient volumes of air into the sample chamber 106 and then disperse that air inside the sample chamber 106, as will be explained shortly. Arrangements other than the motor 118

and fan 116 illustrated in Figure 8B can also be used within the scope of the present invention.

The use of air as the thermal transfer medium, in contrast to other gases and liquids, has the advantages of being inexpensive, readily available, easily mixed, and never making a mess. In the case of the described embodiments, the high surface area-to-volume ratio of the sample containing capillary tubes provides for rapid thermal transfer using air as the thermal transfer medium.

During cooling portions of the thermal cycle, the action of the fan 116 draws ambient temperature air into the housing 102. A vent door 128, articulating on hinge 129, is provided. The vent door 128 is automatically opened by way of a solenoid 132 so that the interior of the housing 102 is sealed off from the upper housing portal 130. In some embodiments, the solenoid 132 is preferably replaced by a stepper motor as is known in the art. The use of a stepper motor allows the vent door 128 to be accurately and incrementally opened and closed in accordance with the needs for heating and cooling the samples. Those skilled in the art will be able to derive an appropriate control mechanism for use with a stepper motor, for example an SC-149 stepper motor controller (available from Alpha Products) as known in the art, using the information set forth herein.

Due to the arrangement of the lower sample chamber portal 120 and the larger cross sectional area and position of the upper sample chamber portal 126, room temperature air is moved into the sample chamber 106 and is dispersed and mixed within the sample chamber 106 by a paddle 124 which is connected to the motor shaft 122. The paddle 124 should rotate at a relatively high rate, for example, fast enough to create air velocities of around preferably about 250, more preferably 500, and most preferably 1000 meters per minute within the sample chamber 106. With the paddle 124, which can be a single or a multivane paddle, rotating at a high speed, air is moved, or drawn, into the sample chamber 106 and vented out of the sample chamber 106 following the path indicated by the dashed line 136. The rotation of the paddle 124 also promotes mixing of the air entering the sample chamber 106 and ensures the most efficient transfer of thermal energy from the surfaces of the thin-walled capillary tubes 108 to the air passing through the sample

chamber 106. It will be appreciated that structures other than those illustrated herein can perform equivalent functions.

As the solenoid 132 is actuated to open the vent door 128, all of the room temperature air moved into the sample chamber 106 is exhausted through a sample chamber upper portal 126 and then through the upper housing portal 130 carrying the heat from the sample chamber 106 to the surrounding atmosphere. The rapid mixing of the air that passes through, and is disbursed in, the sample chamber 106 results in rapid and uniform cooling of the samples.

---

10 Example 2

Figure 9A shows the results of four different temperature/time profiles (A-D) and their resultant amplification products after thirty cycles (A-D). The profiles A and B in Figure 9A were obtained using a prior art heating block device using the prior art microfuge tube. As can be seen in Figure 9A, the transitions between temperatures are slow and many nonspecific bands are present in profiles A and B. Profile B shows improvement in eliminating some of the nonspecific bands (in contrast to profile A) by limiting the time each sample remains at each temperature thus indicating that shorter times produce more desirable results.

Profiles C and D were obtained using the apparatus of Figures 8A-B. As can be seen in Figure 9A, amplification is specific and, desirably, even though yield is maximal in C (60 second elongation) it is still entirely adequate in D (10 seconds elongation).

The optimal times and temperatures for the amplification of a 536 bp fragment of  $\beta$ -globin from human genomic DNA were also determined. Amplification yield and product specificity were optimal when denaturation (93°C) and annealing (55°C) were less than 1 second. No advantage was found to longer denaturation or annealing times. The yield increased with longer elongation times at (77°C) but there was little change with elongation times longer than 10-20 seconds. These unexpected results indicate that the previously available devices used for DNA amplification are not maximizing the conditions needed to optimize the physical and enzymatic requirements of the reaction.

xxx

00531339.030300

Further information can be obtained from: Wittwer, Carl T., Marshall, Bruce C., Reed, Gudrun B., and Cherry, Joshua L., "Rapid Cycle Allele-Specific Amplification with Cystic Fibrosis  $\Delta F_{508}$  Locus," 39 Clinical Chemistry 804 (1993) and Wittwer, Carl T., Reed, Gudrun H., and Ririe, Kirk M., "Rapid DNA Amplification," THE POLYMERASE CHAIN REACTION 174 (1994) which are both now incorporated herein by this reference.

From the information provided in Figure 9A, it can be seen that the embodiments of the present invention subject the samples placed therein to rapid thermal cycling wherein the temperature of the sample is increased and decreased at a rate preferably at least as great as  $0.5^{\circ}\text{C}/\text{second}$ . In the case of the present invention carrying out the polymerase chain reaction, the temperature change is preferably carried out over an approximate range of between  $30^{\circ}\text{C}$  to  $50^{\circ}\text{C}$ . It is preferred that the thermal cycles be carried out quickly enough to complete at least thirty thermal cycles in forty minutes and more preferably complete thirty thermal cycles in twenty minutes and most preferably complete thirty thermal cycles in ten minutes.

The apparatus 100 more preferably increases and decreases the temperature of the sample at a rate at least as great as  $1.0^{\circ}\text{C}/\text{second}$  and even more preferably at a rate at least as great as  $4.0^{\circ}\text{C}/\text{second}$  and most preferably at a rate at least as great as  $10.0^{\circ}\text{C}/\text{second}$ . Critically, the biological sample, not just the surrounding medium and/or the sample container, must undergo the specified thermal change. The previously available devices, while having the drawback of not being able to perform thermal changes as rapidly as the present invention, also did not recognize the problem of changing the temperature of the sample, not just the temperature of the surrounding medium and container, rapidly and uniformly.

Referring now to the chart of Figure 9B, the method of the present invention can desirably achieve thermal cycling preferably at a rate at least as great as  $10^{\circ}\text{C}/\text{sec.}$ , and more preferably at a rate at least as great as  $20^{\circ}\text{C}/\text{sec.}$ , over a temperature range of greater than about  $20^{\circ}\text{C}$ , more preferably over a temperature range of greater than about  $30^{\circ}\text{C}$ , and most preferably over a temperature range of about  $40^{\circ}\text{C}$ . Figure 9B shows the temperature in  $^{\circ}\text{C}$  of the biological sample, not just the surrounding air or container, as the biological sample undergoes thermal cycling. Figure 9B shows a PCR sample beginning at about  $74^{\circ}\text{C}$  and being heated to a

denaturation temperature, indicated at D, of about 92°C for 2 seconds. The sample is then cooled to an annealing temperature, indicated at A, of about 55°C for two seconds. The transition between the denaturation temperature and the annealing temperature covers a range of 37°C in just under 4 seconds providing a rate at least as great as 10°C/sec. The sample is then warmed to an extension temperature of 74°C for five seconds as indicated at E in Figure 9B. The cycling of the sample through the denaturation temperature, the annealing temperature, and the extension temperature is repeated thirty times or as many times as desired.

Figures 9C-G show exemplary cycles of other preferred temperature/time profiles which are achieved by the present invention. It will be understood that those skilled in the art can alter the represented temperature/time profiles to carry out specific processes in accordance with the present invention. Those skilled in the art will also appreciate that the previously available devices and methods, such as devices which conduct heat to and from the sample via a solid or liquid, cannot provide the resulting temperature/time profiles described herein. Moreover, the previously available devices and methods do not suggest or teach the temperature/time profiles described herein. Furthermore, it will be appreciated that the previously available devices and methods utilizing air as the transfer medium, for example previously available chromatographic ovens, cannot provide, and do not suggest or teach, the temperature/time profiles which are described herein and obtained by the practice of the present invention.

In order to provide the fastest thermal cycling time, it is preferred that the lamp (112 in Figures 8A and 8B) be rated at 2000 watts or a plurality of lamps be included which provide similar output. It is also preferred to include a temperature slope control circuit which is represented in Figure 10 in conjunction with an A-bus controller/acquisition system using an 8052 micro controller board with a clock and high level program interpreter available from Alpha Products (model no. SP-127) of Darian, Connecticut. Exemplary programming code used in connection with the described micro controller is included in the Programming Code Appendix A attached hereto and incorporated herein. The programming code provided in Appendix A is a BASIC52 file for serial downloading into the micro controller and provides exemplary temperature slope control during thermal cycling. Use of the 2000 watt heat





5

10

15

20

25

30

5

10

15

20

30

8.3, 2 mM MgCl<sub>2</sub>, 50 µg/ml bovine serum albumin, 0.5 µM each primer, 0.2 mM each dNTP, and 0.4 U native Taq DNA polymerase in a 10 µl reaction. The human beta-globin primers RS42 and KM29 (536 base pairs) are described in C.T. Wittwer, G.C. Fillmore and D.R. Hillyard, "Automated Polymerase Chain Reaction in  
5 Capillary Tubes with Hot Air," Nucl. Acids. Res. 17:4353-4357. Temperature cycling parameters were 94°C for 0 sec., 55°C for 0 sec., and 72°C for 10 sec. Thirty five cycles of amplification were performed with the indicated rates between all temperatures. The samples were electrophoresed on 1.5% agarose gels and stained with 0.5 µg/ml ethidium bromide. Specificity and yield both decrease as the  
10 temperature transition rate decreases.

Fluorescent probes can be used to detect and monitor DNA amplification. As known to those skilled in the art, useful probes include double-stranded-DNA-specific dyes and sequence-specific probes. With the intercalater ethidium bromide, UV-excited red fluorescence increases after amplification. While  
15 microfuge tubes have been used as a sample container for DNA amplification, the embodiments of the present invention described herein advantageously utilize sample containers with many of the characteristics of structures referred to herein as capillary tubes.

The use of the sample containers described herein allows detection of  
20 fluorescence while the sample is held within the container, as will be explained more fully hereinafter. Those skilled in the art will appreciate the number of different schemes of fluorescence detection of DNA amplification which are now available. For example, sequence-specific fluorescence detection is readily possible using the present invention and oligonucleotide hybridization probes. As another example,  
25 dual-labeled fluorescein/rhodamine probes can be cleaved during polymerase extension by 5'-exonuclease activity, separating the fluorophores and increasing the fluorescein/rhodamine fluorescence ratio.

Using the embodiments of the present invention described hereinafter, fluorescence can be measured after temperature cycling is complete, once per cycle as  
30 a monitor of product accumulation, two or more times during a temperature transition, or continuously within each cycle. In contrast to the present invention, previously

09631339-000300

5

15

25

30



Referring now to Figure 11, a schematic view of a preferred rapid temperature cycler with fluorescence detection is provided, generally designated at 300. A forced air hot air source 302 is preferably provided. The forced air hot air source 302 is preferably a commercially available device including a 1600 watt heating coil and fan. A cool forced air cool air source 304 is also preferably provided. The cool forced air source 304 is preferably a 2200 rpm shaded pole blower available in the art from Dayton of Niles, Illinois, model no. 4C006B. It is preferred that the cool air source 304 provide ambient temperature air, but it is within the scope of the present invention to utilize a means for providing fluid that is at a temperature lower than ambient air temperature.

15 In the embodiment of Figure 11, ducts 306 and 308 connect the forced hot air source 302 and the forced cool air source 304, respectively, to a sample chamber 310. The ducts 306 and 308 are preferably corrugated black nylon tubing having a 2.5 cm diameter. The duct 306 is connected to the sample chamber 310 via a port 306A and the duct 308 is connected to the sample chamber 310 via a port 308A.

20 A vent 312 and an exhaust fan 314 function to move air out of the sample chamber 310. Moreover, a means for shielding the interior of the sample chamber 310 from ambient light is integral with the sample chamber 310.

The temperature of the samples within the sample chamber 310 is preferably monitored by a tubular, metal-sheathed thermocouple 316, available from Idaho Technology of Idaho Falls, Idaho, model no. 1844, which is matched in thermal response to the samples held in the preferred sample containers, for example capillary tubes. Importantly, temperature homogeneity within the sample chamber 310 is achieved by mixing the air within the sample chamber 310. It is preferred that such mixing of the air within the sample chamber 310 be carried out by a central sample chamber fan 318. The sample chamber fan preferably includes a 1.7 X 11 cm fan blade available from Idaho Technology, model no. 1862, and a motor available from Idaho Technology, model no. 1861, which creates air velocities of at least 800 to 1000

meters per minute within the sample chamber 310. Such rapid air velocities may not be needed in all applications of the present invention but rapid air velocities promote extensive mixing and temperature homogeneity within the sample chamber 310.

Within the sample chamber 310, a plurality of samples are held in  
5 capillary tubes, some of which are indicted at 320, and are placed in a vertical orientation on a rotatable carousel 322. The carousel 322 is preferably fourteen centimeters in diameter and rotated by a 400 step per revolution stepper motor 324 controlled by a micro stepping drive module 326. The stepper motor 324 is preferably one available from New England Affiliated Technologies of Lawrence,  
10 Massachusetts, model no. 2198364, and the micro stepping drive module 326 is preferably one also available from New England Affiliated Technologies, model no. MDM7 micro stepping drive module, which provides 12,800 steps per rotation of the carousel 322.

Still referring to Figure 11, a fluorescence excitation source 328 is  
15 provided. One preferred arrangement for the excitation path in accordance with the present invention will now be described with one preferred arrangement for the collection path in accordance with the present invention will subsequently be described. The fluorescence excitation source 328 preferably includes a 75 watt xenon arc source 328A focused with an elliptical reflector 328B. The xenon arc  
20 source 328A is preferably available from Photon Technology International of South Brunswick, New Jersey, model no. A1010, with f/2.5 elliptical reflector 328B. The power supply and other components needed to operate the fluorescence excitation source 328 are well known to those skilled in the art. Alternatively, a light emitting diode can be used as a fluorescence excitation source. Those skilled in the art will  
25 appreciate that many different excitation sources can be used within the scope of the present invention.

The radiation emitted by the fluorescence excitation source 328 is focused to about 2 mm using an adjustable iris 334 such as one available in the industry from Rolyn (Covina, California), model no. 75.0125. The light emitted from  
30 the fluorescence excitation source 328 impinges upon a cold mirror 330, which is preferably available from Rolyn, model no. 60.4400, and passes through heat absorbing glass 332, which is preferably one available from Rolyn, model no.

09631339 080300

65.3130. After collimation through a planoconvex lens 336, preferably one available from Rolyn, model no. 10.0260, a 450-490 nm bandpass interference filter 338, preferably one available from Omega Optical of Brattleboro, Vermont, model no. 470RDF40, a focusing planoconvex lens 340, preferably available from Rolyn, model no. 10.0260, and a 1 mm silica window 342, preferably available from Omega, to prevent condensation on the just described optical components during temperature cycling. Using the described excitation path, a 5-7 mm section of one capillary sample tube 320A is illuminated.

Still referring to Figure 11, the collection path for collecting the fluorescence emitted from the sample 320A will be described next. The optics of the collection path include a 1 mm silica window 344 which is placed in the optical path to prevent condensation on the other optical components. Two opposed aspheric lenses 346A&B, preferably available from Rolyn, model no. 17.1175, function to focus emitted fluorescence onto a 2 x 10 mm slit 348. The slit 348 can preferably be fabricated from cutting exposed X-ray film and the slit 348 functions as a spatial filter. After the slit 348 (acting as a spatial filter), the emitted fluorescence is imposed upon a 35 mm electronic shutter 350 operated via an electronic shutter control 352. The 35 mm electronic shutter 350 is preferably a Uniblitz shutter model no. VS35 and the electronic shutter control 352 is preferably driver model no. D122, both available from Vincent Associates of Rochester, New York. A collimating aspheric lens 354, preferably one available from Rolyn model no. 17.1175, is also provided.

A filter 356 is also included when detection of SYBR® Green I emissions is desired. The filter 356 is preferably a 520-580 nm band pass filter, available from Omega as model no. 550RDF60, which is preferably used for single wavelength acquisition. For detection of other emissions, for example, a combination of a dichroic filter 358 and wavelength filters 358A and 358B can be used. For example, for separation of fluorescein and rhodamine emissions, the dichroic filter 358 preferably consists of a 560 nm dichroic filter, preferably available from Omega, model no. 560 DRLP, and a 520-550 nm band pass filter (358A), preferably available from Omega, model no. 535DF30, for detection of fluorescein, and a 580-620 nm band pass filter (358B), preferably available from Omega, model no. 600DF40, for detection of rhodamine. For separation of fluorescein and Cy5 emissions, the dichroic

filter 358 preferably is a 590 nm dichroic filter, available from Omega, model no. 590 DRLP, and filters 358A&B preferably consist of a 520-550 nm band pass filter (358A), available from Omega, model no. 535DF30, for detection of fluorescein, and a 660-680 nm band pass filter (358B), available from Omega, model no. 670DF20, for  
5 Cy5 detection. Those skilled in the art will readily appreciate that the use of other components can be readily implemented using the information set forth herein in order to accommodate other fluorescent wavelengths.

Still referring to Figure 11, after being subjected to the respective filter 358A or 358B, the emitted fluorescence is focused through two planoconvex lenses  
10 360A & 360B, each preferably available from Edmund of Barrington, New Jersey, model no. 32970, and onto photomultiplier tubes 362A and 362B, respectively. The photomultiplier tubes ("PMT") 362A and 362B are preferably available from Hamamatsu of Middlesex, New Jersey, model no. R928, and are each enclosed in a suitable housing including appropriate circuitry, preferably one available from Photon  
15 Technology International, model no. 714, with analog acquisition capabilities. A PMT and data acquisition control module 364 is also preferably provided. Manual PMT shutters 366A and 366B, as known in the art, are also provided.

The forgoing described optical components are preferably five centimeters in diameter and mounted in five centimeter universal lens mounts, such as  
20 those available from Rolyn, model no. 90.0190. As can be carried out by those skilled in the art, many of the necessary structural components were machined from black Delrin™ using techniques known in the industry.

Those skilled in the art will appreciate that the rapid temperature cycler with fluorescence detection 300 can advantageously be constructed using light  
25 emitting diodes (LEDs) and photodiodes in place of similarly functioning components represented in Figure 11. Thus, the function of the fluorescence excitation source 328 can be carried out by light emitting diodes. The photomultiplier tubes 362A&B can also be replaced with photodiodes. Additional information regarding suitable light emitting diodes and photodiodes will be provided later herein. It will be appreciated  
30 that technique sensitivity is limited by background fluorescence, most of which comes from the probes, not the detection system. Significantly, stability is generally more important than absolute sensitivity.



Those versed in the art will appreciate that the rapid temperature cycler with fluorescence detection 300 represented in Figure 11 includes the beneficial characteristics of a fluorimetry device with rapid temperature control, a combination nowhere suggested or taught in the art. PCR can be performed and analyzed during  
5 ten to twenty minutes of temperature cycling. The present invention's combination of 1) fluorescence monitoring within each temperature cycle and 2) analysis of the temperature and time dependence of hybridization provides advantages not otherwise obtainable.

The present invention also makes possible single-color fluorescence  
10 methods to monitor product purity and quantify template during PCR. Dyes that monitor DNA strand status are added to PCR reactions for observation during temperature cycling using embodiments of the present invention.

In order to explain some of the benefits which accrue with the present invention, specific examples using the apparatus represented in Figure 11 will now be  
15 provided. DNA amplification was performed in 50 mM Tris, pH 8.3 (25°C), 3 mM  $\text{MgCl}_2$ , 500  $\mu\text{g/ml}$  bovine serum albumin, 0.5  $\mu\text{M}$  of each primer, 0.2 mM of each deoxynucleoside triphosphate and 0.2 U of Taq polymerase per 5  $\mu\text{l}$  sample unless otherwise stated in the following examples. Also in the following examples, human  
20 genomic DNA (denatured for 1 min by boiling) or purified amplification product was used as DNA template. Purified amplification product was obtained by phenol/chloroform extraction and ethanol precipitation (see D.M. Wallace 1987, Large- and small-scale phenol extractions and precipitation of nucleic acids (as described at p. 33-48, in S.L. Berger and A.R. Kimmel (Eds.), Guide to Molecular Cloning Techniques (Methods in Enzymology, Vol. 152) Academic Press, Orlando),  
25 followed by removal of primers by repeated washing through a Centricon 30 micro concentrator (available from Amicon of Danvers, Massachusetts). Template concentrations were determined by absorbance at 260 nm.  $A_{260}/A_{280}$  ratios of templates were greater than 1.7.

In these examples, primers were synthesized by standard  
30 phosphoramidite chemistry, as known in the art, namely, using Pharmacia Biotech Gene Assembler Plus (Piscataway, New Jersey). The 180 base pair fragment of the hepatitis B surface antigen gene was amplified using primers



10

15

25

30

Instruments, Austin, Texas) and a 12-bit multifunction input/output card 368A (available from National Instruments under the designation AT-MIO-E2) in a PC compatible computer 368 utilizing an Intel® 80486 microprocessor running at a clock speed of 120 MHZ. Analog output channels on the input/output card 368A were used to control the sensitivity, i.e. the PMT voltage, of each of the photomultiplier tubes 362A&B. Analog input channels on the input/output card 368A receive the signals from each of the photomultiplier tubes 362A&B. The PC compatible computer 368, through the input/output card 368A, controls the position, rate and direction of movement of the carousel 322. For example, when multiple capillary sample tubes are loaded, the carousel 322 rapidly positions each capillary sample tube 320 sequentially at a monitoring location (the location represented by capillary sample tube 320A) for a 10 - 100 msec acquisition period. For continuous monitoring of a single capillary sample tube, the capillary sample tube is held in the monitoring position while data is preferably acquired every 200 msec. and is averaged in accordance with well-known techniques. Time, temperature, and preferably two channels of fluorescence are continuously displayed via a monitor 368B associated with the computer 368 as fluorescence vs. cycle number and fluorescence vs. temperature plots.

The carousel 322 should be positioned where maximal fluorescence and signals are acquired. When a single capillary sample tube, such as the capillary sample tube 320A, is monitored the signals are acquired every 200 msec with an integrating time constant set on the photomultiplier tube 362A or 362B, or both, at 50 msec. For multiple sample tubes, the time constant is set at 0.5 msec and the carousel is rotated once to locate the precise position where each capillary sample tube 320 provides the maximum fluorescence in each of the two channels. After positioning the capillary sample tube 320A at a location where maximum fluorescence is obtained, the sensitivity of each PMT 362A&B is adjusted and the carousel rotated again to count and locate the position of all the capillary sample tubes 320 in the carousel 322. When only a signal fluorescence acquisition is desired once each amplification cycle during extension, each capillary sample tube 320 is sequentially positioned on the carousel 322 at the monitoring position for 100 msec. Continuous acquisition for multiple tubes can also be obtained by continuously rotating the

carousel 322. Temperature control programming was based upon, and modified from, a commercial rapid temperature cycler available from Idaho Technology under the trademark Rapidcycler™ using an 8051 cross compiler available from Systronics, Salt Lake City, Utah, designated BCI51 and Dallas development system (also available  
5 from Systronics under the designation DPB2).

In practice, the temperature response of the rapid temperature cycler with fluorescence detection 300 is similar to the response obtained with the embodiment of the present invention disclosed in Figures 8A&B allowing 20-30 second cycles (30 cycles in 10-15 min) as represented in the temperature vs. time  
10 chart of Figure 11A (which shows a few cycles of one preferred temperature profile). When a double strand-specific fluorescent dye is present during amplification, fluorescence generally increases as more double stranded product is made. See R. Higuchi, G. Dollinger, P.S. Walsh, and R. Griffith, 1992, "Simultaneous Amplification and Detection of Specific DNA Sequences," Bio/Technology 10:413-  
15 417.

Moreover, it will also be appreciated that double strand specific dyes such as ethidium bromide or SYBR® Green I can be used as generic indicators of amplification. SYBR® Green I dye is preferred over ethidium bromide in many applications because it has an excitation maximum near fluorescein and often  
20 provides a stronger signal with DNA than visible excitation of ethidium bromide.

Fluorescence also depends on temperature, a confounding effect during temperature cycling that is usually eliminated by considering fluorescence once per cycle at a constant extension temperature. However, if temperature, time, and fluorescence are acquired every 200 msec during rapid cycle amplification, a three  
25 dimensional spiral is shown on the monitor 368B as represented in Figure 12. The three dimensional plot represented in Figure 12 is also projected in Figure 12A as a two dimensional plot of temperature vs. time, projected in Figure 12B as a two dimensional plot of fluorescence vs. time, and projected in Figure 12C as fluorescence vs. temperature. The temperature vs. time projection of Figure 12A repeats each  
30 cycle and provides essentially the same information as set forth in Figure 11A. Because fluorescence varies inversely with temperature, the fluorescence vs. time projection shown in Figure 12B at early cycles is a scaled mirror image of the

00631339-080300



5

10

25

30

Figure 14, and subsequent figures, for different initial template copy number. When the data are normalized as the percent maximal fluorescence of each capillary sample tube 320, one hundred initial copies are clearly separated from ten copies. However, the difference between one and ten copies is marginal, and no difference is observed between zero and one average copies per capillary sample tube 320.

Double strand dyes depend on the specificity inherent in the amplification primers. As will be appreciated by those skilled in the art, nonspecific amplification at high cycle numbers can limit detection sensitivity to about one hundred initial template copies (see Figure 14). With rapid cycling taught by the present invention, further improvements in amplification specificity are obtained further improving the overall DNA amplification performance.

Quantitification with sequence-specific probes has a similar dynamic range as double stranded DNA dyes but, as shown in the plots of Figures 15A and 15B, appear to discriminate even a single initial template copy from negative controls.

When low copy number detection and quantification are needed, additional specificity is provided by fluorescent probes that require hybridization for signal generation. Cleavage of a dual-labeled exonuclease probe is one technique which is capable of distinguishing a single template copy from a negative control as shown by Figure 15. Figure 15 show plots of fluorescence ratio vs. cycle number for different initial template copy number, according to the legend provided in Figure 14A.

Signal generation with 5'-exonuclease probes is dependent not only on DNA synthesis, but requires hybridization and hydrolysis between the fluorophores of the dual-labeled probe. This hydrolysis reduces quenching and the fluorescence ratio of fluorescein to rhodamine emission increases. For more information on this technique, see L.G. Lee, C.R. Connell and W. Bloch, 1993, "Allelic Discrimination by Nick-translation PCR with Fluorogenic Probes," Nucl. Acids Res. 21:3761-3766 & Livak, K.J., S.J.A. Flood, J. Marmaro, W. Giusti and K. Deetz, 1995, "Oligonucleotides with Fluorescent Dyes at Opposite Ends Provide a Quenched Probe System Useful for Detecting PCR Product and Nucleic Acid Hybridization," PCR Meth. Appl. 4:357-362).

Figure 25 shows fluorescence PCR results from a probe with five intervening bases between fluorescein and rhodamine labels. The forty-five cycle amplification was completed in 20 minutes using the rapid temperature cycler with fluorescence detection 300 of Figure 11. By monitoring the fluorescence ratio once  
5 per cycle, a  $10^9$  fold range of initial template concentration could be distinguished. The amplification curves are shifted approximately 3-4 cycles for each 10-fold change in initial template concentration.

Although the final fluorescence signal is decreased when low copy numbers are amplified (presumably because of decreased amplification efficiency),  
10 quantification between zero and one hundred copies is readily possible. The signal generated by exonuclease probes is cumulative and only indirectly related to product concentration. Hence, the fluorescence signal continues to increase even after the amount of product has reached a plateau. Using the information contained herein, those skilled in the art can formulate appropriate standards to control for efficiency of  
15 amplification and cleavage in order to carry out absolute quantification.

Fluorescence vs. temperature plots of 5'-exonuclease probes confirm that probe hydrolysis is the mechanism of signal generation. In Figure 16, a fluorescence vs. temperature plot of two-temperature cycling is shown with the  $\beta$ -actin exonuclease probe. In each cycle the fluorescence ratio varies linearly with  
20 temperature and there is little, if any, hysteresis. The signal increases each cycle during the annealing/extension phase when probe hydrolysis occurs. Although the fluorescence of both fluorescein and rhodamine decreases with increasing temperature (data not shown in the figures), the rate of change is greater for rhodamine, resulting in an increasing ratio with increasing temperature. No temperature-dependent  
25 hybridization effects are apparent with the 5'-exonuclease probe.

In contrast, when the fluorescence signal is dependent only on hybridization, fluorescence ratio vs. temperature plots show a different pattern with hysteresis during two-temperature cycling, as plotted in Figure 17. The plots in Figure 17 represent the results obtained using two adjacent hybridization probes which are  
30 present, an upstream probe labeled 3' with fluorescein and a downstream probe labeled 5' with Cy5™. The probes are separated by a 1 base pair gap. During the annealing/extension phase of the reaction, the probes hybridize resulting in



fluorescence from resonance energy transfer. The temperature dependence of hybridization can be used to detect mutations by a shift in the melting curve.

From the foregoing discussion, it will be appreciated that fluorescence monitoring during DNA amplification is an extraordinarily powerful analytical technique. Using the rapid temperature cycler with fluorescence detection 300, productive and cost efficient real time monitoring, sequence-specific detection, and quantification can be achieved in five to twenty minutes, depending on the number of initial template copies present.

15 Furthermore, the system and results represented in Figures 11-17 is particularly suited for continuous monitoring of a biological reaction using fluorescent dyes. For example, with precise temperature control and double-strand-specific dyes, product purity can be estimated by melting curves. With rapid temperature control provided by the present invention, absolute product concentration can be determined  
20 by reannealing kinetics. The present invention advantageously provides rapid temperature changes and strict intra-sample temperature homogeneity which is not available in the prior art. In contrast to the prior art, the present invention utilizes sample containers with a high surface area to volume ratio, (for example by using the preferred capillary sample tubes 320 in Figure 11) and uses air as the thermal transfer  
25 medium providing rapid control of sample temperature not otherwise obtainable. For example, sample temperature vs. time plots obtained when processing samples in the sample containers of the present invention show sharp spikes at denaturation and annealing temperatures (showing rapid temperature response) in contrast to the prior art conical plastic tubes which require several seconds for all of the sample to reach  
30 thermal equilibrium. Moreover, the sample containers of the present invention provide improved results over using etched silicon or glass chips as sample containers since the thermal cycle times and thermal homogeneity of the present invention are



5

20

30

15 In accordance with the present invention and as is possible using the apparatus of the present invention, rapid temperature cycling with minimal annealing and denaturation times improves quantitative PCR and increases the discrimination of allele specific amplification. Rapid cycling for cycle sequencing reduces sequencing ambiguities and minimizes "shadow banding" in dinucleotide repeat amplifications.

20 In accordance with the present invention, for long PCR up to 35 kb, yield is improved when the sample is exposed as little as possible to high denaturation temperatures.

In contrast to the previous approach to PCR which treat PCR as three reactions, denaturation, annealing, extension, each of which occur at three different temperatures (as represented in Figure 18A), one aspect of the present invention provides that a kinetic paradigm for PCR renders important improvements. Using a kinetic paradigm for PCR (as represented in Figure 18B), the temperature vs. time curve consists of continuous transitions between overlapping reactions. The method and apparatus of the present invention is particularly efficient at carrying out PCR under the kinetic paradigm. Figure 18C is a graph representing different time/temperature profiles near an annealing temperature of 55°C. In Figure 18C, the solid trace shows a centrally positioned “spike” representing the temperature of response of a 10 µl sample. In contrast, the traces shown as short and long line

segments in Figure 18C represent the temperature responses of samples obtained using heat block instruments. As can be seen from Figure 18C, the embodiments of the present invention produce annealing segment "spikes," with the advantages discussed herein, in contrast to the temperatures "plateaus" according to the conventional wisdom in the art.

The previously available instrumentation used for detection presented many drawbacks. Rapid, precise temperature cycling is provided by the system of the present invention described herein, in contrast to previously available instrumentation that is five to ten times slower. With the continuous fluorescence monitoring also provided by the system of the present invention, the temperature dependence of hybridization can be followed. By following hybridization during temperature cycling, the number of probes and/or spectral colors required can be minimized. That is, different products and mutations can be detected by their dynamic melting characteristics, rather than going to the trouble of synthesizing different fluorophore-labeled probes for each DNA species that is to be detected.

In order to provide an embodiment of the present invention that is most cost effective, a high intensity light emitting diode is used instead of a xenon arc source or a laser for sample illumination, and photodiodes are used for detection. Samples are loaded into glass capillary sample tubes, or alternatively into composite glass/plastic sample containers (see Figure 21A-D) in a 96-well format that does not require heat sealing. The present invention thus provides real time fluorescence analysis in a cost effective manner. Real time fluorescence control of temperature cycling improves amplification quality. For example, if the temperature of samples is increased only until denaturation occurs, product exposure to high temperatures is minimized. This increases yield by limiting product and enzyme degradation and increases specificity by limiting amplification of products with a high melting temperature.

Reference will next be made to Figure 19, which provides a diagrammatic representation of another preferred embodiment of the present invention configured for continuous monitoring of a single sample. It will be understood, however, that the structures represented in Figures 19 and 20 can also be incorporated into a system which automatically processes multiple samples, such as the apparatus

represented in figure 11 and as will be explained shortly herein. In the embodiment of Figure 19, a single sample holder 402 is placed in a holding bracket 404 positioned at the intersection of a temperature-controlled air stream and a linear optical path. The sample holder 402 includes a tube 402A which has many of the desirable characteristics of a capillary tube. In accordance with the present invention, different configurations of capillary tubes can be used and the tube 402A preferably has a rectangular cross section. The biological sample preferably is held at a bottom end of the tube 402A as indicated at 402B. A cap 402C is also preferably provided on the sample holder 402.

Reference will next be made to Figures 19A-19E which compare the effect of different configurations of sample containers on the temperature response of the sample itself. The temperature-time tracings shown in Figure 19E correspond to the response obtained using the sample container configurations represented in Figures 19A-C, respectively. Figure 19D represents a sample container which is less preferred for use in the present invention and is included for comparison. Using the information set forth herein, those skilled in the art can arrive at optimum sample container configurations for particular applications of the present invention. Further information regarding each of the sample container configurations represented in Figures 19A-D are set forth below.

<u>Figure</u>	<u>Surface Area</u> (mm <sup>2</sup> /10μl)	<u>Fluid Column</u> <u>Length</u> (mm)	<u>Sample Volume</u>	<u>Source</u>
19A	77	47	10μl	Kimble KIMAX #46485-1
19B	42	13.8	34μl	Kimble KIMAX #46485-15
19C	32	8	59μl	Kimble KIMAX #34500-99

<u>Figure</u>	<u>Surface Area</u> (mm <sup>2</sup> /10μl)	<u>Fluid</u> <u>Column</u> <u>Length</u> (mm)	<u>Sample</u> <u>Volume</u>	<u>Source</u>
19D	18	N/A	10μl	MICROAMP <sup>™</sup> tube of Perkin-Elmer Cetus GeneAmp PCR System 9600

In the apparatus of Figure 19, an excitation radiation source 418, preferably an LED and most preferably a blue LED, is provided to illuminate the sample holder 402. The radiation emitted by the excitation radiation source 418 passes through aspheric focusing lenses 420 and an excitation bandpass filter 422 and the radiation is focused onto the sample holder 402.

The optical components illustrated in Figure 19 are preferably held in an optical housing 412. A housing 406 is also provided. A fan 408 is provided to move air through an air duct 414 and over the sample holder 402 held in the sample bracket 404. A temperature unit 410 is placed in the air flow path to provide heating or heating and cooling for the air passing over the sample holder 404. A nozzle 416 effectively directs the air over the sample holder 404.

The emissions which are given off by the sample pass through two more aspheric lenses 420 and an emission bandpass filter 424 and are received by a photo detector 426, which preferably is a photo diode. Those skilled in the art can readily provide the control components needed to advantageously operate the apparatus represented in Figure 19 using the information set forth herein.

Figures 19F and 19G are side and end views, respectively, of one preferred sample container 403 which utilizes a rectangular capillary tube 403A. The capillary tube 403A is preferably one available from Vitro Dynamics Inc. having dimensions of 1mm X 3mm X 50mm. A first cap member 403B and a second cap member 403C are held together by a screw 403D, the screw 403D also functioning as a holder for the capillary tube 403A.

Figures 19H and 19I, respectively, show two possible orientations of a rectangular capillary tube 403A when detecting fluorescence of the sample contained therein. Figure 19H shows the rectangular capillary tube 403A oriented so that its edges are in line with the optical axis of the excitation and detection optics ("edge excitation and detection"). Figure 19I shows the rectangular capillary tube 403A oriented so that its faces are in line with the optical axis of the excitation and detection optics ("face excitation and detection"). Surprisingly, the fluorescence signal obtained from the edge detection orientation shown in Figure 19H is about three-fold to about five-fold higher than obtained with the face detection orientation shown in Figure 19I. The desirable characteristics of using the edge detection orientation shown in Figure 19H is at least partially due to total internal reflection which takes place in the capillary tube 403A which concentrates the fluorescence signal to the extremities of the capillary tube 403A.

Figure 20 shows the optical components of another preferred embodiment in accordance with another aspect of the present invention. The optical components represented in Figure 20 are preferably incorporated into the thermal cycling and sample handling structures represented in Figure 21, which will be more fully described shortly, but which can also be used with many different arrangements to provide monitoring (most preferably continuous monitoring) of a sample undergoing the polymerase chain reaction.

In contrast to the arrangements previously disclosed herein, the optical excitation and detection paths are combined in the embodiment of Figures 20 and 21, referred to herein as an epifluorescent path, rather than a linear path. In the embodiment of Figures 20 and 21, the excitation and emission radiation follow the same optical path between the capillary tube and the dichroic element used in the excitation path. A capillary tube is particularly adapted for use in the embodiment of Figures 20 and 21 due to the total internal reflection (also referred to as "light piping") along the length of the capillary sample tube which is exploited to increase both excitation and emission intensities.

In the embodiment of Figures 20 and 21, to accommodate maximal light piping, the optical axis is parallel to the length of the capillary tube (paraxial) with the tip of the capillary tube positioned at the focal point. Assuming a refractive

index of about 1.33 for the sample being detected, about 12.3% of emitted light is guided to the tip. It is understood that centrifuge action can be used to move the sample to the tip of the capillary tube.

Figure 22A charts the effectiveness of light piping when detecting fluorescence at the tip of the capillary tube and shows a 10-fold increase in signal intensity by viewing the tip (closed diamonds) rather than the side (open circles) of the capillary sample container. Also, as indicated in Figure 22B, the results obtained using capillary sample tubes of two different sizes and which were filled to different lengths with dsDNA stained with SYBR® Green I are plotted. As can be surmised from Figures 22A and 22B, the observed epifluorescence increases as more sample is added to the tube, although the fluorescence efficiency decreases.

The optical properties of the emission from a capillary were investigated by stimulating fluorescence in a capillary filled with a fluorescein solution at 470 nm. The emission from a blunt end of the capillary was seen to be homogenous across the face of the capillary as opposed to concentrated in the glass as would be expected if the emission were the result of evanescent wave fluorescence.

The optical components represented in Figure 20 carry out paraxial epifluorescent illumination of the capillary tip, which provides advantageous results not otherwise obtainable. In Figure 20, an excitation radiation source 468 is preferably a blue LED, such as one known in the industry as a super bright LED and available from LEDtronics. The emitted fluorescence signals are acquired by photo detectors 466A and 466B. The excitation radiation source 468 and the photo detectors 466A and 466B are supported on a mounting board 468 which also includes necessary circuitry and which integrates filters with the photo detectors 466A and 466B. A preferred mounting board is available from Ealing Electrooptics which includes 0.5 inch interference filters with high performance silicon photodiodes in TO5 packages. The excitation and detection components are supported directly on the mounting board 468 with associated electronics. It is preferred that the optical components are preferably  $\leq 1.0$  inches in diameter. A collimating lens 454, two dichroic filters 456A and 456B, a mirror 458, interference filters 460A-C, and aspheric focusing lenses 462A-C direct the radiation to and from the sample.

While the embodiment of the present invention represented in Figure 20 utilizes only two colors/wavelengths when performing an analysis, those skilled in the art can readily adapt the embodiment to provide three, or more, color analysis. To provide three or more color analysis, the apparatus represented in Figure 20 can accommodate additional dichroic filters and photo detectors. Moreover, it is within the scope of the present invention to allow simultaneous separation of wavelengths onto a linear photo detector array, as is available in the industry, for multicolor acquisition. When a linear photo detector array is used in accordance with the present invention, it is preferred that a prism or diffraction grating be utilized in cooperation with a lens and a photo detector array or CCD for detection of multiple wavelengths. One preferred linear photo detector array available in the industry collects 15-30 wavelength bins of 10-20 nm each between 500 and 800 nm. Various configurations of optical components, for example the Littrow autocollimating configuration for gratings used in most monochrometers, can be arrived at using the information set forth herein to arrive at the best accommodation between collection efficiency, spectral resolution and spatial requirements. The apparatus of Figure 20 will now be further described incorporated into an automated thermal cycling apparatus represented in Figure 21.

Figure 21 provides a schematic representation of another presently preferred embodiment 400 of the present invention which includes rapid temperature cycling components, sample handling components, and the optical components represented in Figure 20, all working together to provide fluorescence detection at the tip of the sample containers (epifluorescence). The rapid temperature cycler with epifluorescence detection 400 represented in Figure 21 provides particular advantages. It is to be understood that this described embodiment is merely exemplary of the present invention and that those skilled in the art can arrive at many different arrangements for carrying out the invention claimed herein.

In the embodiment represented in Figure 21, air is taken in through an aperture 470 and generally follows the flow path indicated by the lines 472. The temperature of the air, and thus the temperature of the plastic/glass sample container 450, is preferably adjusted using a 400 watt heating cartridge 474 which is preferably one available from Reheat, Inc. The heating cartridge 474 is positioned within a



central duct 476. A fan 498 is provided to move the air in the indicated path 472. The fan is driven via a shaft 496 and a motor 494. The motor 494 is preferably a DC rare earth brush motor which is preferably available from Escap AG. and having a maximum rpm of 15,000. When heating the plastic/glass sample tubes 450, the heating cartridge is proportionally controlled and the fan is run at a relatively low speed (12 volts, 0.5 amp) to provide temperature homogeneity for all of the plastic/glass sample containers 450. When cooling the plastic/glass sample containers 450, the heating cartridge 474 is disabled and the motor 494 is run at a fast speed (for example with the above-mentioned preferred motor maximum speed is obtained by applying 27 volts, 1.4 amps). The fan 498 forces air into the aperture 470 and out via exhaust ports 471.

In the preferred rapid temperature cyler with epifluorescence detection 400, it is preferred that twenty-four plastic/glass sample containers 450 (two of which are represented in Figure 21) be symmetrically arranged around the heating cartridge 474 and the central duct 476. The plastic/glass sample containers 450 are received by sleeves 451 which (due to their offset structure) allow for precise adjusting of the position of the individual plastic/glass sample containers 450 in a circular carousel 480. The sleeves 451 are preferably fabricated from brass. The off-axis structure of the sleeve 451 allows each sleeve 451 to be aligned so that the tip of the glass/plastic sample container 450 can be precisely adjusted to be at the optical focal point represented in Figure 21, both laterally and longitudinally, at the time that the rapid temperature cyler with epifluorescence detection 400 is fabricated.

The carousel 480 is supported on a bearing 482 above a housing 490. The carousel 480 is positioned by a stepper motor 488 provided with a drive gear 484 connected to the motor 488 via a shaft 486. The stepper motor 488 is microstepped (using a controller (not explicitly represented in Figure 21) from New England Affiliated Technologies) to provide over 10,000 steps per revolution of the carousel 480, providing precise positioning of each the plastic/glass sample containers 450. The interior of the housing 490 is provided with an insulative material 492, preferably in accordance with the previously described insulative material. Baffles 476 function to form the exhaust port 471 and to block ambient light.

Figures 21A-D provide additional detailed views of the plastic/glass sample containers 450 and will be referred to for an explanation of the preferred method of using the same. The plastic/glass sample container 450 includes a capillary tube portion 450B which is closed at one end. The capillary tube portion 450B can take many different configurations and is not limited to only a capillary tube type structure. It is, however, preferred that the volume of fluid held by the plastic/glass sample containers 450 be not more than 1 milliliter in order to promote sample temperature homogeneity and rapid thermal cycling. For example, it is preferred that the material from which the capillary tube portion 450B is fabricated have a thermal conductivity in the range from about 20 to about 35 in accordance with the formula.

$$\frac{\text{cal cm}}{\text{cm}^2 \text{ s degree C}} \times 1$$

Further information regarding the thermal conductivity of different glasses can be obtained from R.C. Weast & M.J. Astle, HANDBOOK OF CHEMISTRY AND PHYSICS, page E-6 (1982)(CRC Press) which is now incorporated herein by reference. The plastic/glass sample containers 450 are also provided with a reservoir portion 450C which is preferably fabricated from an appropriate plastic and joined to the open end of the capillary tube portion 450B. While many different materials can be used for the reservoir portion 450C, it is preferred that a plastic material be formed in a funnel-like shape and attached to the capillary tube portion 450B.

20 A sample S is loaded into the composite plastic/glass sample container  
450 using a pipette P, or some other appropriate instrument, either manually or using  
an automated process. It is preferred that the volume of the sample be in the range  
from about .01  $\mu\text{l}$  to about 10,000  $\mu\text{l}$ , more preferably in the range from about .01  $\mu\text{l}$   
to about 100  $\mu\text{l}$ , and most preferably in the range from about .01  $\mu\text{l}$  to about 10  $\mu\text{l}$ ,  
25 with about 5  $\mu\text{l}$  being the most preferred volume. Once a sample has been added to  
each plastic/glass sample container 450, the plastic/glass sample containers 450 are  
centrifuged at low speed to place the samples at the tips of the closed end of the  
capillary portion 450B, so that the sample forms a 0.2 - 2.0 cm column of fluid 450A  
as represented best in Figure 21B. A stopper 450D (which is preferably configured as  
30 a plastic plug) is then placed in the reservoir portion 450C to seal the plastic/glass

sample container 450 as shown best in Figure 21C and the plastic/glass sample container 450 is placed in the sleeve 451 in the rapid temperature cyclers with epifluorescence detection 400. It is also within the scope of the present invention to provide different structures to seal the capillary tube portion 450B.

5           The capillary tube portion 450B of the glass/plastic sample container 450 is preferably a glass capillary tube available in the industry having 0.8 mm inner diameter and a 1.0 mm outer diameter, and which is closed/sealed on one end and flared at the other end for receiving the plastic reservoir 450C. The glass/plastic sample containers 450 can be readily and economically fabricated. The shape of the  
10 tip 450E of the capillary tube portion 450B is optimized for optical efficiency. Flat tips as well as tips with various outside curvatures and inside curvature are all contemplated within the scope of the present invention. Those skilled in the art can select the most efficient configuration for the tip.

          As can be discerned from Figures 21A-D, the addition of plastic  
15 loading and sealing structures to a capillary tube provides great advantages and allows efficient use of glass capillary tubes while retaining their desirable thermal characteristics. It will be appreciated that it is within the scope of the present invention to add the samples to the plastic/glass sample containers 450, and to subject the samples to centrifuging, in a 96-well format. Moreover, it is within the scope of  
20 the present invention to load the plastic/glass sample containers individually into the rapid temperature cyclers with epifluorescence detection 400 and it is also within the scope of the present invention to provide an embodiment of the present invention to load the plastic/glass sample containers 450 in a 96-well format or some other format.

          Advantageously, the composite plastic/glass sample containers 450  
25 provide a convenient, inexpensive sample holder. With the embodiment of Figure 21, it is preferred that fluorescence is acquired from single samples one to ten times each second. When acquiring fluorescence from multiple samples at the preferred rate, the samples need to be moved into position by rotation of the carousel 480 relatively rapidly. With the preferred stepper motor 488 and appropriate control devices (which  
30 can be selected using the information contained herein) each of the twenty-four samples can be rapidly and accurately moved to the monitoring position represented in Figure 21.

When the fluorescent signal from each sample is acquired for 100 msec., the signal variation (with repositioning) is <1%. It will be appreciated that it is within the scope of the present invention to decrease the signal acquisition time, increase the transit speeds, and also observe the coefficient of variation from repeated  
5 sampling. When twenty-four samples are processed, and the carousel is rotated without stopping at a rate between one and ten revolutions per second, each sample has 0.37-3.7 msec of excitation and detection.

Using the information set forth herein, one skilled the art can select whether the fluorescent signal is integrated via software or hardware. In one preferred  
10 embodiment, a graphical programming language is used in connection with the rapid temperature cycler with epifluorescence detection 400, such as one known in the industry as LabView (available from National Instruments), which has subprograms for peak detection and integration. In another preferred embodiment, integration is done in hardware with variable integration time (user adjustable sensitivity control) so  
15 that the signals reach a level optimal for analog-to-digital conversion.

Using the rapid temperature cycler with epifluorescence detection 400 represented in Figure 21, continuous monitoring of the sample as the reaction is ongoing allows determination of temperature cycling requirements during amplification, based on continuous observation of annealing, extension, and  
20 denaturation. This is in contrast to the prior art where all cycling parameters are determined and programmed before amplification begins. In accordance with the prior art, using complementary oligonucleotides equivalent to the lowest melting primer, the annealing efficiency is controlled even during early cycles. In many cases, extension and denaturation can only be monitored with dsDNA dyes during later  
25 cycles when enough product has been made. Significantly, such a requirement is not usually a problem because denaturation and extension conditions are made permissive enough to amplify most products, and data from the first amplification can be used to optimize subsequent runs.

Still referring to Figure 21, a user interface and instrument control 500  
30 can be fabricated using the information set forth herein in connection with the embodiment of Figure 11. As one preferred example of a user interface and instrument control 500, a PENTIUM™ microcomputer running the LabView

programming language with a 12-bit multifunction input/output card (available from National Instruments) provides data acquisition and control. It is preferred that the analog output signals be used to adjust the amplifiers associated with the photo detectors 466A and 466B. Analog input channels also measure the temperature of the samples via a thermocouple 499 as well as the fluorescent detected from the sample by the photodiodes. The user interface and instrument control 500 represented in Figure 21 also provides digital I/O control of the excitation radiation source 468, the direction of the stepper motor 488, the heating cartridge 474, and the fan 498.

When continuous fluorescence monitoring of PCR samples containing the dsDNA dye SYBR Green I or fluorescently labeled oligonucleotide probes can be used to monitor hybridization and melting during individual amplification cycles. This information can be used by preferred arrangements for the user interface and instrument control 500 to provide improved and customized thermal cycling conditions. The benefits of using hybridization information for temperature cycling include:

(A) Ensuring that complete denaturation of the PCR product occurs with each cycle while:

Minimizing exposure to excessively high denaturation temperatures, thus, avoiding heat induced damage to the amplification products and polymerase.

Increasing reaction specificity by minimizing the denaturation temperature which selects against products with a  $T_m$  higher than the intended amplification product.

(B) Maximizing the amplification efficiency by ensuring adequate time for product extension with each cycle while:

Minimizing the amount of time required for amplification by allowing no longer than needed to complete product extension.

Enhancing reaction specificity by selecting against products longer than the intended amplification product.

(C) Maximizing the amplification efficiency by ensuring adequate time for product extension each cycle while:

Minimizing the amount of time required for amplification by allowing no longer than needed to complete product extension.

Enhancing reaction specificity by selecting against products longer than the intended amplification product. These would require longer than the allotted time to complete product extension.

(D) Initiating thermal cycling changes dependent on the level of fluorescence obtained or the current efficiency of amplification. For example, over-amplification and nonspecific reaction products can be minimized by terminating thermal cycling when the efficiency drops to a certain level. As another example, temperature cycling can be modified to initiate slower temperature ramps for melting curve acquisition when the fluorescence becomes detectable. This saves time because the slower ramps need not be used on earlier cycles. Other desirable changes may become evident on continued practice of the invention.

(E) Minimizing over-amplification damage to PCR product and/or initiation of melting curve acquisition before over-amplification has increased the background of nonspecific reaction products.

In accordance with the present invention, the user interface and instrument control 500 can follow preprogrammed time/temperature set points and/or, advantageously, can acquire detected fluorescence values and then use the acquired detected fluorescence values to alter or adjust one or more reaction parameters in real time to optimize the results obtained. As used herein, the term “reaction parameter” includes, but is not limited to, any parameter which is used as a basis for controlling a reaction. Such reaction parameters include, but are not limited to, denaturation temperature and time, primer annealing temperature and time, probe annealing temperature and time, enzyme extension temperature and time, and number of cycles. In general, control of the reaction is initially based on an estimate of reaction parameters from the fluorescence data. The original fluorescence data is either acquired as a change in fluorescence over time (temperature specific rates of denaturation, annealing, and extension), a change in fluorescence over temperature (product or probe  $T_m$ ), or a change in extent of amplification (amplification yield and efficiency). These rates,  $T_m$ 's and their first and second derivatives are used to

5 between those carried out by a portion of the user interface and instrument control 500 (which preferably can be an IBM compatible computer using programing based upon the teachings set forth herein)(Blocks 500A-500E in Figure 22C) and those carried out by the remaining components (Blocks 500A, and 500G-500S in Figure 22C) of the rapid temperature cyler with epifluorescence detection 400. It is to be understood  
10 that the block diagram of Figure 22C is merely exemplary and many different arrangements can be used to carry out the present invention.

sample containing the reaction mixture at the temperature at which the product is seen to have completely melted. Each cycle of the reaction uses this fluorescence value as a target. The approach being described in this example uses two stages in to provide a time lag to accommodate the requirement of sending the fluorescence values to a separate PC computer. With each product melting step, the temperature is increased until the fluorescence reaches an intermediate value, then the power applied to the heating device is reduced so that a temperature ramp of approximately 3 °C per second is imposed so that the PC computer has adequate time to analyze the fluorescence and convey to other components that product denaturation has occurred. The resulting time/temperature plot is shown in Figure 22D. Figure 22D shows a characteristic increase in the melting temperature after twenty cycles as the concentration of amplification product grows. This is due to the fact that product  $T_m$  is a function of product concentration.

30 at a combined annealing/extension temperature, the fluorescence of the sample is  
monitored and this information is used to ensure that adequate, but not excessive, time  
had been allowed for product extension. Fluorescence is monitored at ten second

Figure 22E shows a fluorescence/time plot which exhibits a characteristic increase in the dwell time at the combined annealing/extension temperature as the concentration of amplification product grows. This is due to the fact that as the primer concentration and polymerase become limiting more time is needed to complete product extension with each cycle.

As a yet another example of the advantages of the arrangement shown in Figure 22C, amplification plateau will be discussed. At the end of each annealing/extension step, the fluorescence value is acquired and stored. When this value increases to 1.2 times the lowest end-cycle fluorescence value and had subsequently stopped increasing below a user settable ratio (typically 1.00 - 1.02) the thermal cycle is terminated. Alternatively, a melting curve accusation step is initiated by entering a slow 0.1 °C to 0.2 °C/second temperature ramp through the product  $T_m$  and monitoring the fluorescence of the sample continuously. The resulting fluorescence/time plot shown in Figure 22D shows that after twenty-five cycles of amplification the ratio of cycle-by-cycle fluorescence growth fell below 1.00 and the reaction terminated. It will be appreciated that this approach can be used to acquire a high resolution melting curve for each sample. As a sample reaches its amplification plateau, a melting curve can be acquired for that sample, then regular temperature cycling can resume until another reaction reaches its amplification plateau.

Figure 22E illustrates useful temperature vs. time segments for  
25 fluorescence hybridization monitoring. Product melting curves are obtained during a  
slow temperature increase to denaturation. By quickly lowering the temperature after  
denaturation to a constant temperature, product, probe, or primer annealing can be  
detected. Probe melting curves are obtained by slowly heating through temperatures  
around the probe  $T_m$ . Those skilled in the art can readily utilize the system  
30 represented in Figure 21 to provide the necessary analysis, in real time if desired,  
during temperature cycling to provide heretofore unavailable information on the



characteristics of the product, probe, and primer using the hardware and software described herein.

Absolute quantification of product is also advantageously carried out in accordance with the present invention. Continuous monitoring of double stranded DNA formation allows direct, absolute DNA quantification by reannealing kinetics. The sample temperature is quickly dropped from the denaturation temperature and held constant at a lower temperature that is still high enough to prevent primer annealing. The rate of product reannealing then follows second order kinetics. When different concentrations of DNA are tested, the shape of the reannealing curve is characteristic of the DNA concentration (see Figure 26). For any given PCR product and temperature, a second order rate constant can be measured. Once the rate constant is known, any unknown DNA concentration can be determined from experimental reannealing data. The curves can be fit by non-linear least squares regression during temperature cycling in real time using the LabView programming environment (explained previously). Cooling is not instantaneous, and some reannealing occurs before a constant temperature is reached, but regression analysis allow for this in accordance with the present invention. (see Figure 27). The technique requires pure PCR product, but this can be verified by melting curves also obtained during temperature cycling. Quantification by reannealing kinetics is independent of signal level and not affected by sample volume differences.

Figure 28 is a schematic representation of another embodiment of the present invention which includes many of the structures included in the embodiment of Figure 21. In order to provide a succinct description of the embodiment of Figure 28, only those significant differences between those components represented in Figure 21 and those components represented in Figure 28 will be explained with the understanding that one skilled in the art can readily use the information contained herein to fabricate embodiments of the present invention. Figures 27A and 27B are cross sectional schematic views of the embodiment represented in Figure 28 in a run mode and a load mode, respectively.

The embodiment of Figure 28 is a rapid temperature cycler, generally designed at 502, with fluorescence detection at the tip of the sample containers with automatic positioning of the sample containers in two dimensions which improves the

fluoresce signal which is obtained from the sample. Figure 29 is a perspective view of the exterior of the embodiment of the present invention including the components illustrated in the schematic representation of Figure 28.

As seen in both Figures 28 and 29, a removable circular sample tray  
5 483 holds thirty-two samples. The removable circular sample tray 483 is placed into the rapid temperature cyclor 502 so that it engages a carousel 481 which is driven by a motor 488. As the carousel 481 rotates, a hall effect position locator is used to precisely position the carousel 481 so that the each sample is precisely positioned over a flourimeter assembly 459. The flourimeter assembly 459 preferably includes a LED  
10 source 459A, three photodiodes 459B, focusing lenses 459C, and a filter assembly 459D. The flourimeter assembly 459 is similar in structure and function to that represented in Figure 20.

Most advantageously, the fluorimeter is mounted on a slicer bearing  
493 which is moved by a lateral stepper motor 491. As the carousel 481 rotates, the  
15 composite plastic/glass sample containers 450 are precisely positioned over the fluorimeter assembly 459 in the direction of the carousel and the position is noted by the apparatus via the hall effect position locator 495 while the lateral stepper motor 491 adjusts the position of the fluorimeter assembly 459 is adjusted in a second dimension, and the position noted. Thus, the rapid temperature cyclor 502 provides  
20 for improved placement of a plurality of samples into the apparatus using a removable sample tray 483 and provides for improved detection of a fluorescence signal from a sample.

Provided in Figures 30A-V are detailed schematic diagrams showing the preferred configuration of the electrical components of the rapid temperature  
25 cyclor 502 represented in Figures 28 and 29. It is to be understood that the diagrams of Figures 30A-V are merely one preferred arrangement for carrying out particular aspects of the present invention and these diagrams are not intended to be limiting of the scope of the present invention. In order to improve the clarity of the diagrams, the notations which are commonly used in the industry are maintained on these diagrams  
30 and are referenced in the corresponding parts list provided below.

Parts List - MAIN				
	<u>Item</u>	<u>Quan-tity</u>	<u>Reference</u>	<u>Part</u>
5	1	1	BT1	3V LITHIUM
	2	9	C1,C2,C3,C8,C9,C13,C18,C24,C26	.1UF
	3	7	C4,C5,C10,C12,C14,C15,C17	1UF
	4	2	C7,C6	100UF
	5	6	C11,C16,C19,C20,C21,C22	18pF
10	6	1	C23	47UF
	7	2	C25,C27	22UF
	8	2	C28,C29	10UF
	9	1	F1	1A
	10	1	IC1	AD594
15	11	2	IC2,IC3	DS5000FP
	12	1	IC4	LM324
	13	8	IC5,IC7,IC10,R13,R17,R18,R21,R22	10K
	14	2	IC6,IC8	MS62256
	15	2	IC9,IC10	DS2003
20	16	1	IC11	TLC1451
	17	1	IC12	7432
	18	1	IC13	PT5101
	19	1	IC14	PT5102
	20	1	IC15	7404
25	21	1	IC16	PIC16C54
	22	1	IC17	MAX232
	23	1	IC18	LM4040
	24	1	IC19	LTC1293
	25	1	IC20	LTC1286
	26	1	IC21	LM385 1.2
	27	1	IC22	LTC1144

Parts List - MAIN			
Item	Quan-tity	Reference	Part
28	2	IC23,IC24	PVG612S
29	1	JP1	HALL SENSOR
30	1	JP2	FLO1
31	1	JP3	FLO3
32	1	JP4	FLO2
33	1	JP5	MAIN HEADER
34	1	J1	CON2
35	1	LED1	STEP
36	1	LED2	TEMP
37	2	LED3,LED4	RED/GREEN LED
38	1	P1	SERIAL CONNECTOR
39	1	Q1	2N5484
40	10	Q2,Q3,Q4,Q5,Q6,Q7,Q8,Q9,Q10,Q11	NDS351
41	1	R1	4.87K 1%
42	4	R2,R4,R5,R6	10K 1%
43	1	R3	2.74K 1%
44	1	R7	200
45	8	R8,R9,R10,R11,R19,R20,R28,R29	470
46	2	R15,R12	100
47	3	R14,R16,R23	1K
48	4	R24,R25,R26,R27	4.7K
49	1	S1	TYPE J
50	1	Y1	20.0000
51	2	Y3,Y2	14.745600

## 5

10

15

20

25

Parts List - POWER BOARD

<u>Item</u>	<u>Quan-tity</u>	<u>Reference</u>	<u>Part</u>
28	1	R3	360
29	7	R11,R13,R14,R16,R17,R18,R19	10K
30	1	R12	4.7K
31	1	R15	1K
32	1	R20	261
33	1	R21	866
34	1	R22	650
35	1	R23	180
36	2	S1,S2	110/220
37	1	T1	TRANSFORMER FLAT COMPACT
38	1	VR1	LM2575

Parts List - INTEGRATOR

<u>Item</u>	<u>Quan-tity</u>	<u>Reference</u>	<u>Part</u>
1	2	C1,C2	1UF
2	1	C3	.01
3	1	IC1	ACF2101
4	1	IC2	OPT301
5	1	IC3	OPA627
6	1	IC4	REF200
7	1	J1	CON6
8	1	P1	500
9	1	R1	30M
10	1	R2	100K
11	2	R3,R4	10K

5

20

In accordance with another embodiment of the present invention a handling system is provided for loading small volume sample vessels with liquid samples, particularly samples to be analyzed by detection of emitted fluorescence. The sample vessel typically has a volume of less than 1 ml, and it can be in the form of a tube (i.e. a capillary tube) or a "flat capillary" wherein the capillary space is defined by two spaced-apart plates or sheets sealed along their edges. The sample vessel typically has a volume to external surface area ratio of about 1 mm<sup>3</sup>, more typically less than about 0.5 mm<sup>3</sup>. Capillary tubes having an inner diameter of less than 1 mm have a volume to surface area ratio of less than 0.25 mm<sup>3</sup>. The vessel used in accordance with the present invention is preferably formed from an optically transparent material. Preferred materials are optically transmissible for light having a wavelength ranging from about 400 to about 800 nm. The use of such material will allow the detection of a fluorescent signal generated in a liquid sample held by the vessel. Moreover the use of vessels with a low volume to surface area ratio for analyzing fluorescence from a

fluorescent sample enables more efficient detection of the fluorescence due to enhanced total internal reflection.

Vessels having a high surface area to volume ratio (or conversely, a low volume to surface area ratio) can be difficult to load with liquid samples.

- 5 Advantageously, the sample handling system of the present invention helps to overcome such difficulties. In accordance with one embodiment a vessel having a high surface area to volume ratio and an open end is provided with a funnel cap that fits onto the open end of the vessel to facilitate loading of liquid samples into the vessel. The funnel cap includes a first sample receiving port and a second sample transfer port and means for releasably fixing the funnel cap on the vessel so that the sample transfer port of the funnel cap and the open end of the vessel are in alignment. In one embodiment the funnel cap is of plastic or rubber construction and is formed so that the inner diameter of the sample transfer port frictionally engages the outer diameter of the vessel proximal to its open end. However, other means of coupling the funnel cap to the vessel are known to those skilled in the art and are within the scope of the invention, including the use of adhesives, clamps, clasps and the like. In one embodiment the sample handling system further comprises a plug for frictional fit sealing engagement with the sample receiving port of the funnel cap. However any device or material that effectively seals the opening of the funnel to prevent contamination or evaporation of the loaded sample is suitable for use with the present invention.

- Advantageously the vessels of the present invention can be used in a method for enhancing detection and efficiency of acquisition of fluorescence in a sample comprising a fluorophore. The method comprises the steps of placing a sample in a vessel having walls composed of an optically transparent material and defining a volume having at least first and second dimensions. The first dimension is less than the second dimension and the ratio of volume to external surface area of the vessel is less than 1mm. Enhanced detection and efficiency of acquisition of fluorescence generated from the sample is achieved by detecting fluorescence along an axis substantially parallel to a wall along the second dimension of the vessel. In one embodiment, sample fluorescence is induced by fluorophore-excitatory illumination of the sample wherein the sample is illuminated along an axis



5

10

20

30

which are indicated at 9. The sample vessel port 7 and sample passageway 8 are formed for receiving and fixing sample vessel 9 to the disc 2. In one embodiment the sample vessel 9 is releasably fixed to the carousel 1 to allow the removal of the sample vessel and its replacement with another sample vessel to allow for multiple use of the carousel 1. In an alternative embodiment the sample vessels 9 are permanently fixed to, or formed as an integral component of, the disc 2. In one embodiment the sample vessel 9 is fixed to the disc 2 by frictional contact between the sample vessel 9 and at least a portion of the sample passageway 8 proximal to said sample vessel port 7. Other conventional means for fixing the sample vessel in communication with the sample vessel can be used. For example, complementary screw threads can be formed on the surface of the sample passageway 8 and on the exterior surface of the sample vessel 9. In addition adhesives or any other fixation means known to those skilled in the art can be used in accordance with the present invention to fix the sample vessel 9 to the disc 2. The top and bottom surfaces of the carousel of the present invention are preferably formed to allow multiple carousels to be stacked one on top of another so that a stack of multiple carousels can be releasably engaged with a motor drive shaft and rotated simultaneously as a unit as shown in Figure 32.

The embodiment shown in Figure 32 includes a stepper motor 504 and a drive shaft 506 which functions to hold and rotate the carousels generally indicated at 1. A chamber fan 508 is used to generate the air flow indicated by the arrows 512. A heating device 510 functions to heat the air which passes by the sample vessels 9. A fluorimeter assembly 514 includes an LED source 514A, photodiodes 514B, focusing lenses 514C, and a filter assembly 514D. A fluorimeter stepper motor 516 functions to move the fluorimeter assembly 514 in the direction of arrow 518. Those skilled in the art can readily fabricate embodiments of the present invention fashioned after the arrangement represented in Figure 32 using the information set forth herein.

In another embodiment (not shown) the carousel comprises a disc having a top surface, a bottom surface, an outer edge extending therebetween, a sample receiving port in the top surface, a sample vessel port in the bottom surface and a sample passageway communicating with said sample receiving port and the sample vessel port. The sample vessel port and sample passageway are formed for

In one embodiment the sample passageway of the disc comprises a first portion having a central axis substantively parallel to the top and bottom surfaces of the disc and a second portion having a central axis forming an acute angle with the top and bottom surfaces of the disc. In this embodiment the sample vessel port and sample passageway are formed for receiving and fixing a sample vessel to the disc such that the sample vessel extends from the disc at an acute angle relative to the bottom surface of the disc.

Carousel 1 is further provided with means for closing the sample receiving ports 6A, 6B, and 6C. The closure means can be a plug (not shown) that fits into the sample receiving port 6 and frictionally engages the adjacent walls of the sample passageway, or for example, adhesive backed tape, for application to the top surface to effectively seal the opening of the sample receiving port to prevent contamination or evaporation of a loaded sample. Carousel 1 is releasably engaged with a drive shaft for rotation. Any suitable engagement means well known to those of ordinary skill in the art can be used including frictional engagement, or the use of screws, bolts, locking pins or clamps. In one embodiment, the disc 2 is formed as ring having a center hole formed for receiving a drive shaft (see 506 in Figure 32). The end of the drive shaft is preferably provided with structures for holding the discs 2 in place.

The carousel 1 of the present invention can be used to deliver a liquid sample to a sample vessel 9. In one embodiment the sample vessel 9 is a capillary vessel containing a predetermined mixture (for example a reagent mixture) that interacts with one or more components of the introduced sample. In accordance with one embodiment the predetermined mixture is added to the sample vessel before positioning a capillary sample vessel into the sample vessel port. Alternatively the sample vessel is prepackaged with a predetermined mixture. The predetermined mixture may comprise reagents that react or interact with the sample to produce a detectable signal or to produce a derivative product.

The sample passageway 8 of the carousel 1 are optionally provided with one or more barriers 10 that prevent a liquid sample delivered through sample

receiving ports 6A, 6B, and 6C from flowing to the sample vessel port 7 absent a biasing force on said liquid sample. The term "barrier" is used herein to include any structure that impedes the free flow of a liquid sample delivered into a sample receiving port to the sample vessel port. Examples of suitable barriers for use in the sample passageway of the carousel of the present invention include depressions or wells formed in the sample passageway, sample passageway narrowing projections or annular rims that extend from the surface of the sample passageway, porous membranes, directional valves, or flaps that are biased in a closed position.

The barriers are formed so that the liquid sample can overcome the barrier by application of a biasing force on a liquid sample present in the sample passageway and blocked by the barrier. The application of biasing force on the sample is preferably provided by the centripetal force generated by rotation of the carousel. Therefore, in a carousel having a plurality of sets of sample receiving ports 6A, 6B, and 6C in the top surface, each set with a corresponding sample passageway and sample vessel port, samples can be added individually to the various sample receiving ports and the barrier will localize the liquid sample and prevent the samples from flowing to the respective sample vessel ports. After all of the samples are delivered into the respective receiving ports, the carousel is rotated to deliver the samples to the respective sample vessel port and into an attached sample vessel.

In accordance with one embodiment, each sample passageway of the carousel communicates with a single sample vessel port and a plurality of sample receiving ports. In accordance with that embodiment, the sample passageway can optionally include a central passageway that branches to communicate with multiple sample receiving ports, or alternatively, as illustrated in Figures 31A&B multiple sample receiving ports 6A, 6B, and 6C are aligned along a common axis that extends radially from the center of the disc, each of said ports communicating through one passageway with a sample vessel received in the sample vessel port. The sample passageway can be provided with one or more barriers 9A that prevent a sample added to any one of the plurality of sample receiving ports from flowing to the sample vessel port absent a biasing force on said liquid sample. Furthermore, each sample passageway can be provided with multiple barriers, each of which require a different amount of biasing force to transfer a sample over the barrier. In accordance with this

embodiment, after delivery of the samples to the respective sample receiving ports, individual samples can be selectively transferred to the sample vessel port and into the sample vessel by controlling the rate of rotation of the carousel.

For example, a first sample can be delivered into a first sample receiving port and a second sample can be delivered to a second sample receiving port wherein the first and second sample receiving ports communicate with a common passageway and the first and second sample receiving ports are each provided with a barrier that prevents flow of the respective first and second sample. The barriers allow the disc to be provided as part of a kit with predetermined amounts of selected reagents, catalysts, enzymes, oils, etc. being preloaded into the sample passageway via one or more of the sample receiving ports.

In one embodiment the barrier for the second sample receiving port is formed so that a greater biasing force must be applied to the sample delivered to the second sample receiving port to pass its associated barrier than is required for a sample delivered to the first sample receiving port to pass its associated barrier. In accordance with this embodiment, rotation of the carousel at a first rate will deliver the first sample to the sample vessel port and into the sample vessel, while the second sample is prevented from flowing to the sample vessel port and into the sample vessel. Rotation at a increased second rate will then enhance the centripetal force on the second sample and result in the delivery of the second sample to the sample vessel port and into the sample vessel. Based on this principle, different samples can be delivered to multiple sample vessel ports that communicate with a common passageway and after all the samples have been loaded, the individual samples can be delivered to the sample vessel port and into the sample vessel one at a time or simultaneously by controlling the rate of rotation of the carousel. In one embodiment a first sample, comprising a fluorophore is added to a first sample vessel port and a second sample comprising oil is delivered to the second vessel port. The carousel is rotated to deliver the first sample into the sample vessel followed by the oil. The oil (or another liquid that effectively seals the first sample within the sample vessel) functions both to decrease evaporation of the first sample and to reduce the risk of contamination of the first sample.

One preferred sample vessel carousel of this invention includes three sample receiving ports preferably, but optionally, arranged in radial alignment and in fluid communication with a common sample vessel. In accordance with this embodiment, about 1 to about 5  $\mu\text{l}$  of an oil overlay, preferably dyed black, is present in prepackaged form, or delivered to the radially innermost sample receiving port. The oil overlay comprises mineral oil and about 0.01% to about 1% organic black dye such as Waxoline® Black OBP available from Zenica, Inc. of Wilmington, DE. About 1 to about 9  $\mu\text{l}$  of a reagent master mix is present in prepackaged form or is delivered to the radially outer most sample receiving port. The reagent master mix comprises a portion of, or all the necessary reaction components. A liquid sample containing the template nucleic acid to be tested is delivered manually or robotically into the radially intermediate sample receiving port. The disc is then rotated at a rate

One preferred sample vessel carousel of this invention includes three sample receiving ports preferably, but optionally, arranged in radial alignment and in fluid communication with a common sample vessel. In accordance with this embodiment, about 1 to about 5  $\mu\text{l}$  of an oil overlay, preferably dyed black, is present in prepackaged form, or delivered to the radially innermost sample receiving port. The oil overlay comprises mineral oil and about 0.01% to about 1% organic black dye such as Waxoline® Black OBP available from Zenica, Inc. of Wilmington, DE. About 1 to about 9  $\mu\text{l}$  of a reagent master mix is present in prepackaged form or is delivered to the radially outer most sample receiving port. The reagent master mix comprises a portion of, or all the necessary reaction components. A liquid sample containing the template nucleic acid to be tested is delivered manually or robotically into the radially intermediate sample receiving port. The disc is then rotated at a rate

that transfers the sample to the reagent compartment, but at a rotated rate insufficient to deliver the mixture into the sample vessel. The sample and reagent can optionally be mixed by rapid changes in the rate of the rotation of the disc. The disc is then rotated at a higher rate that causes the sample and reagent mixture, but not the oil, to move into the sample vessel. The disc is then rotated at still a higher rotation rate to deliver the oil overlay to the sample vessel. The oil will overlay the aqueous sample because of its lower density and will block light passage because of its dye content. The selective transfer of oil, reagents and sample by altering the rate of carousel rotation is achieved by a combination of: 1) varying the diameter of the fluid communication passageways; 2) varying the size or shape of the physical barriers present in the fluid communication passageways; and 3) by using the dependence of centrifugal force on the varying distance (radius) of each sample receiving port from the center of the disc.

The carousel of the present invention can be releasably engaged with the drive shaft and a motor (506 and 504, respectively in Figure 32) for rotating the carousel. Furthermore, individual carousels of this invention can be stacked upon one another and engaged with a drive shaft for simultaneous rotation (as shown in Figure 32). In accordance with another aspect of the present invention a device is provided for monitoring the fluorescence of a sample held within a sample vessel (see 514 in Figure 32). The sample vessel comprises an optically transparent material and has walls defining a volume having at least first and second dimensions wherein the first dimension is less than the second dimension and wherein the ratio of volume to external surface area of the vessel is less than 1mm. In one embodiment the device comprises a chamber, a sample vessel holder, a light emitting source mounted in said chamber and positioned to illuminate the sample vessel along an axis substantially parallel to a wall along the second dimension of the vessel and a light detector mounted in said chamber and positioned to measure fluorescence from the sample vessel along an axis substantially parallel to a wall along the second dimension of the vessel. The light emitting source and the light detector in accordance with one embodiment are mounted on a platform that can be raised and lowered (as indicated by arrow 518 in Figure 32). In this embodiment, the light emitting source and the light detector can be positioned to measure fluorescence from the sample vessels



10

15

25

30



5  
10  
15  
20  
25  
30

25           The carousel of the present invention can also be used for delivering a liquid sample to a capillary sample vessel. The carousel comprises a disc having a top surface, a bottom surface and an outer edge extending therebetween, a sample receiving port in the top surface, a sample vessel port in the outer edge and a sample passageway communicating with the sample receiving port and the sample vessel

30   port. The sample vessel port and the sample passageway are formed for receiving and fixing a sample vessel to the disc. The method of using the carousel to deliver a liquid sample to a capillary sample vessel comprises the steps of selecting a carousel

The present invention is also directed to a system for detecting the presence of a target nucleic acid sequence in a sample. The system comprises a pair of oligonucleotide probes that hybridize to adjacent regions of the target nucleic acid sequence, wherein one of said probes is labeled with an acceptor fluorophore and the other probe labeled with a donor fluorophore of a fluorescence energy transfer pair. Preferably, the donor fluorophore emission and the acceptor fluorophore absorption overlap less than 25%, the acceptor fluorophore has a peak extinction coefficient greater than  $100,000 \text{ M}^{-1}\text{cm}^{-1}$  and upon hybridization of the two probes with the target sequence, the donor and acceptor fluorophores are within 15 nucleotides of one another. In another embodiment the donor fluorophore emission and the acceptor fluorophore absorption overlap less than 20% and upon hybridization of the two probes with the target sequence, the donor and acceptor fluorophores are within 5 nucleotides of one another, and more preferably within 3 nucleotides of one another.

In view of the foregoing, it will be appreciated that the present invention provides an apparatus for accurately submitting biological samples to thermal cycling and for quickly and accurately varying the temperature of biological samples, most advantageously adjusting one or more reaction parameters in real time or according to a predetermined temperature versus time profile. The present invention also provides an apparatus suitable for subjecting a number of different biological samples to rapid thermal cycling and also provides a thermal cycling apparatus having a thermal transfer medium of low thermal mass which can effectively subject samples to a large temperature gradient over a very short period of time.

Moreover, the present invention provides an apparatus which can subject a biological sample to rapid thermal cycling using air as a thermal transfer medium and which provides a system and method for performing PCR rapidly and for simultaneously monitoring the reaction. Still further, the present invention also provides a system and method for performing PCR rapidly and also continuously

5 January 31, 1995 which is now incorporated herein in its entirety.

The present invention may be embodied in other specific forms without departing from its spirit or essential characteristics. The described embodiments are to be considered in all respects only as illustrative and not restrictive. The scope of the invention is, therefore, indicated by the appended claims rather than by the foregoing description. All changes which come within the meaning and range of equivalency of the claims are to be embraced within their scope.

PATENT APPLICATION

of

Carl T. Wittwer,

Kirk M. Ririe,

Randy P. Rasmussen

and

David R. Hillyard

for

## CONTAINER FOR CARRYING OUT AND MONITORING BIOLOGICAL PROCESSES

## PROGRAMMING CODE APPENDIX A

```

1  NEW      GOTO 50:REM      ***I/O scanner (I=*(0), Q=*(72))***
2  10      IF (KEY(55307).AND.64)=0 THEN( Z(0)=KEY(55306) : GOTO 3
3  3        IF : WHILE (KEY(55307).AND.128)=128
4  4          I=0 : DO : I=I-1 : KEY(55306)=ASC$(0,I)
5  5          WHILE ASC$(0,I)<>11 : RETURN
6  6          DO : WHILE (KEY(55307).AND.64)=64
7  7          I=0 : DO : Z(I)=KEY(55306)
8  8          I=I-1 : UNTIL Z(I-1)=32 : RETURN
9  9
10 10      REM      ***error check***
11 15      IF TX<0 THEN TX=0
12 16      RETURN
13 20      IF TX<-151 THEN TX=-151
14 21      RETURN
15 25      IF TX>151 THEN TX=151
16 26      RETURN
17
18 30      GOSUB 50:TF1=TF:TI=TIME
19 32      DO:UNTIL TIME-TI>0.2
20 34      GOSUB 50:TF2=TF:T2=TIME
21 36      SL=(TF1-TF1)/(T2-T1):IF SL<>0 THEN IT=(TG-TF2)/SL
22 38      PRINT USING(###.##), 'SL',SL, 'IT=',IT, 'AT TEMP=',TF2
23 40      IF IT>1 GOTO 30
24 42      RETURN
25
26 49      REM      ***temperature measurement***
27 50      B=KEY(55304):F=KEY(55303):G=KEY(55304):TF=(F*256+G+100)/38.3: RETURN
28
29
30 60      CLEAR : STRING 150,72 : DIM Z(72)
31 72      DTEMP=35:DTEMP1=DS=10
32 74      ATEMP=60:ATEMP1=AS=20
33 76      EXTEND=74:EXTEND1=ES=1:CYCLE=100
34 80      PRINT USING(###), '*',TAB (20), 'TEMP',TAB (30), 'TIME',TAB (40), 'SLOPE'
35 82      PRINT '1) DENATURATION', TAB(30), DTEMP,TAB(30), DTEMP1,TAB(40), DS
36 84      PRINT '2) ANNEALING', TAB(30), ATEMP,TAB(30), ATEMP1,TAB(40), AS
37 86      PRINT '3) EXTENSION', TAB(30), EXTEND,TAB(30), EXTEND1,TAB(40), ES
38 88      PRINT '4) CYCLES', CYCLE:PRINT
39 110     INPUT 'TYPE # TO CHANGE OR 0 IF OK (CHOICES ARE 0,1,2,3, OR 4)?',J
40 114     IF J=0 THEN GOTO 170
41 116     ON J GOSUB ,110,110,150,150
42 118     PRINT:GOTO 80
43 120     INPUT 'ENTER DENATURATION Temp, Time, Slope: ',DTEMP,DTEMP1,DS
44 125     IF DTEMP>38 GOTO 110
45 126     RETURN
46 130     INPUT 'ENTER ANNEALING Temp, Time, Slope: ',ATEMP,ATEMP1,AS
47 140     IF AS<0 GOTO 110
48 142     IF AS>11 GOTO 110
49 145     RETURN
50 150     INPUT 'ENTER EXTENSION Temp, Time, Slope: ',EXTEND,EXTEND1,ES:RETURN
51 160     INPUT 'ENTER Number of Cycles:',CYCLE:RETURN
52 170     S(0)='start 100;sp1 100;hfl 1;tl1 1;ext': GOSUB 3 : GOSUB 7
53 175     FOR I=0 TO 10 : PRINT CHR(Z(I)), : NEXT I : PRINT
54
55 200     S(0)='sol md:hfl 2;wfo =x.:sl':GOSUB 3: REM open
56 210     PRINT 'WAIT FOR COOLING':KEY(55304)=3
57 220     DO : GOSUB 30 : UNTIL TP<45
58 230     PRINT 'TEMP=',TP, ' LOAD SAMPLES AND THEN PRESS A KEY'
59 240     W=GET : IF W=0 GOTO 240
60
61 300     TMC=0 : W6=0: CLOCK 0 : TIME=0 : DBY(71)=0 : CLOCK 1
62 305     S(0)='hfl 1;sp1 100;sol;wfo =x...:sl':GOSUB 3: REM CLOSE
63 310     FOR V=1 TO CYCLE
64 315     PRINT 'D', : TG=DTEMP:KEY(55309)=4:IF V=1 THEN TG=DTEMP+(DS/4)

```





# THE 1950S

Carl T. Wittwer,

Randy P. Rasmussen

David R. Hillyard

for

## CONTAINER FOR CARRYING OUT AND MONITORING BIOLOGICAL PROCESSES

## PROGRAMMING CODE APPENDIX B



Temperature processor(IC 2) hex file (Temp24a.hex):

:0300000020033C8  
:05000B00C0D002312508  
:05001300C0D00230D74F  
:05002300C0D00232BA5A  
:05002B00C0D0023042CC  
:20003300C2AFC2D4C2D3752100752200752300D215D28AAF81758130858115C0D0C0F0C098  
:20005300E0C08753C6FB1230897527405388CF5389F012317D78FC12305D12324912365439  
:20007300C20EC20FC208759850D2ACD2AFE587D0E0D000D001D0D030E20530E4028012206A  
:20009300E602800DBF070280089038A974FA02373975C7AA75C755438740907471C082C03F  
:2000B3008378187A00C000C0029037A6C082C083D003D001D005D004D002D00012309590D6  
:2000D3007084C082C08378007A00D083D082E8F0907060C082C08378007A00D083D082E8C1  
:2000F300F0A3EAF002013BD281D29090705CC2D57400C4C0A8C2AFC282C2907908C282A28C  
:2001130081339280D282D9F5F5F07904C282A28133D282D9F7C282D290D0A8F0E5F0C454C9  
:200133000FA3F030D501222202015BD2A3C2A122D2A3D2A122C2A3D2A122C2A222D2A22248  
:20015300D2A022C2A02222212015612014812015012015A907055C082C08378017A00D0FF  
:2001730083D082E8F078007A001230D278007A0012319D123216A8887A00C000C00278F7F2  
:200193007A00C000C002D003D001D002D00012322D8888A8887A00C000C00278407A00C06B  
:2001B30000C002D003D001D002D000123231888890232FA882AA83123235D217907471127D  
:2001D30033931233D8123416D297907087C082C08378007A00D083D082E8F0122328122357  
:2001F3002890705AC082C08390705CE0F8A3E0FAD083D082E8F0A3EAF01225FF907087E09D  
:20021300F87A00C000C00278007A00C000C002D003D001D002D00012342FE84A6003020110  
:20023300F112232812015312014390707EC082C08378007A00D083D082E8F0A3EAF09070B2  
:200253007CC082C08378007A00D083D082E8F0A3EAF0907088C082C08378017A00D083D07B  
:2002730082E8F078007A001230D2907074C082C08378017A00D083D082E8F0A3EAF0907025  
:2002930062C082C083907060E0F8A3E0FAD083D082E8F0A3EAF0907060C082C08390708947  
:2002B300C082C083907074E0F8A3E0FAC000C00278087A00C000C002D003D001D002D00099  
:2002D30012348EC000C002C000C00278E87A03C000C002D003D001D002D000123508E84A0D  
:2002F30070030230C278017A00C000C002D003D001D002D0001234CFC000C002D003D0018E  
:20031300D002D00012350EC000C002D003D001D002D0001234CFC000C002D003D0018E  
:200333007080C082C083907089C082C083907074E0F8A3E0FAC000C00278087A00C000C062  
:2003530002D003D001D002D0001234CFC000C00278017A00C000C002D003D001D002D000F0  
:2003730012350EC000C002C000C00278E87A03C000C002D003D001D002D000123508E84AEB  
:2003930070030230C278017A00C000C002D003D001D002D0001234CFC000C002D003D001ED  
:2003B300D002D00012350EC000C002D003D00189828B83123528C000C00278D27A00C0007F  
:2003D300C002D003D001D002D0001234CFC000C002907089C082C083907074E0F8A3E0FA94  
:2003F300C000C00278087A00C000C002D003D001D002D0001234CFC000C00278027A00C05B  
:2004130000C002D003D001D002D00012350EC000C002C000C00278E87A03C000C002D00336  
:20043300D001D002D000123508E84A70030230C278017A00C000C002D003D001D002D00093  
:200453001234CFC000C002D003D001D002D00012350EC000C002D003D00189828B831235D1  
:2004730028C000C002D003D001D002D00012350ED083D082E8F0A3EAF0907085C082C08320  
:20049300907089C082C083907074E0F8A3E0FAC000C00278087A00C000C002D003D001D000  
:2004B30002D00012348EC000C00278037A00C000C002D003D001D002D000123512C000C06B  
:2004D30002C000C00278E87A03C000C002D003D001D002D000123508E84A70030230C27880  
:2004F300017A00C000C002D003D001D002D0001234CFC000C002D003D001D002D000123582  
:200513000EC000C002D003D00189828B83123528D083D082E8F0907089C082C0839070740D  
:20053300E0F8A3E0FAC000C00278087A00C000C002D003D001D002D00012348EC000C002B9  
:2005530078057A00C000C002D003D001D002D000123512C000C002C000C00278E87A03C0CF  
:2005730000C002D003D001D002D000123508E84A70030230C278017A00C000C002D003D060  
:2005930001D002D0001234CFC000C002D003D001D002D00012350EC000C002D003D00189C4  
:2005B300828B83123528C000C00278C87A00C000C002D003D001D002D00012352DE84A700F  
:2005D300030205DA020A16907068C082C083907088E0F87A00C000C00278017A00C000C046  
:2005F30002D003D001D002D000123512D083D082E8F0A3EAF0907068C082C083907068E018  
:20061300F8A3E0FAC000C002907089C082C083907074E0F8A3E0FAC000C00278087A00C0BD  
:2006330000C002D003D001D002D00012348EC000C00278047A00C000C002D003D001D0025B  
:20065300D000123512C000C002C000C00278E87A03C000C002D003D001D002D00012350866

00000000"00000000

**DEBILITATION**



:2014B300342FE84A70030214BF120143907072C082C08378017A00D083D082E8F0A3EAF002  
:2014D300907085C082C08378007A00D083D082E8F0907054C082C08378007A00D083D08210  
:2014F300E8F090706AC082C08378017A00D083D082E8F0A3EAF0907080E0F8A3E0FAC00090  
:20151300C00278007A00C000C002D003D001D002D00012354DE84A70030215FD78007A00FD  
:201533001230D278007A0012319D12232812189A12359CC000C002907080E0F8A3E0FAC097  
:2015530000C002D003D001D002D000123508C000C002907087E0F87A00C000C00278017A51  
:2015730000C000C002D003D001D002D000123596C000C002D003D001D002D00012356EC076  
:2015930000C00290707EE0F8A3E0FAC000C00278007A00C000C002D003D001D002D0001255  
:2015B3003596C000C002D003D001D002D00012356EC000C00290707CE0F8A3E0FAC000C0FD  
:2015D3000278007A00C000C002D003D001D002D000123596C000C002D003D001D002D00097  
:2015F30012356EE84A600302153D907057E0F87A00C000C00278007A00C000C002D003D0F8  
:2016130001D002D00012352DC000C00290707EE0F8A3E0FAC000C00278007A00C000C00255  
:20163300D003D001D002D000123596C000C002D003D001D002D00012356EC000C002907075  
:2016530087E0F87A00C000C00278017A00C000C002D003D001D002D00012342FC000C0026A  
:20167300D003D001D002D00012356EE84A7003021863907057E0F87A00C000C00278027A1B  
:2016930000C000C002D003D001D002D00012342FC000C002907057E0F87A00C000C00278D5  
:2016B300037A00C000C002D003D001D002D00012342FC000C002907056E0F87A00C000C0B3  
:2016D3000278017A00C000C002D003D001D002D00012342FC000C002D003D001D002D000FD  
:2016F30012356EC000C002D003D001D002D000123575E84A70030217C0907057C082C08344  
:2017130078017A00D083D082E8F012232812189A907087E0F87A00C000C00278017A00C017  
:2017330000C002D003D001D002D00012342FC000C00290707EE0F8A3E0FAC000C0027800CA  
:201753007A00C000C002D003D001D002D000123596C000C002D003D001D002D00012356EDA  
:20177300C000C00290707CE0F8A3E0FAC000C00278007A00C000C002D003D001D002D000C7  
:20179300123596C000C002D003D001D002D00012356EE84A600302171D90707CC082C08310  
:2017B30078007A00D083D082E8F0A3EAF0907057C082C08378017A00D083D082E8F01225A7  
:2017D300F412232812189A907087E0F87A00C000C00278017A00C000C002D003D001D0029B  
:2017F300D00012342FC000C00290707EE0F8A3E0FAC000C00278007A00C000C002D003D0A3  
:2018130001D002D000123596C000C002D003D001D002D00012356EC000C00290707CE0F8E2  
:20183300A3E0FAC000C00278007A00C000C002D003D001D002D000123596C000C002D003AA  
:20185300D001D002D00012356EE84A60030217D490707EC082C08378007A00D083D082E849  
:20187300F0A3EAF090707CC082C08378007A00D083D082E8F0A3EAF078007A001230D2787D  
:20189300007A0012319D22907057E0F87A00C000C00278037A00C000C002D003D001D002A1  
:2018B300D00012342FE84A7003021993907068C082C083907062E0F8A3E0FAC000C00290C7  
:2018D300705CE0F8A3E0FAC000C002D003D001D002D000123537D083D082E8F0A3EAF09004  
:2018F3007068C082C083907068E0F8A3E0FAC000C00290706AE0F8A3E0FAC000C002D00325  
:20191300D001D002D00012348ED083D082E8F0A3EAF0907068E0F8A3E0FAC000C0027801BB  
:201933007A00C000C002D003D001D002D000123553C000C002907056E0F87A00C000C0020B  
:2019530078007A00C000C002D003D001D002D00012342FC000C002D003D001D002D000126B  
:20197300356EE84A7003021990907056C082C08378017A00D083D082E8F01225F4021CBA13  
:20199300907085E0F87A00C000C00278007A00C000C002D003D001D002D00012342FC000EC  
:2019B300C002907085E0F87A00C000C00278C87A00C000C002D003D001D002D00012342F02  
:2019D300C000C002D003D001D002D000123575E84A7003021A04907062C082C08390706064  
:2019F300E0F8A3E0FAD083D082E8F0A3EAF0021CBA907058C082C08312359CC000C00290DB  
:201A13007085E0F87A00C000C002D003D001D002D0001234CFC000C00278207A00C000C07B  
:201A330002D003D001D002D0001234CFD083D082E8F0A3EAF0907058C082C083907058E027  
:201A5300F8A3E0FAC000C002780A7A00C000C002D003D001D002D0001231AAD083D082E83E  
:201A7300F0A3EAF0907068C082C083907058E0F8A3E0FAD083D082E8F0A3EAF0907058C03A  
:201A930082C08312318FC000C00278647A00C000C002D003D001D002D00012350EC000C027  
:201AB3000278197A00C000C002D003D001D002D0001231AAC000C002907085E0F87A00C038  
:201AD30000C002D003D001D002D0001234CFC000C00278087A00C000C002D003D001D00262  
:201AF300D0001234CFC000C00278647A00C000C002D003D001D002D0001231AAD083D082BC  
:201B1300E8F0A3EAF0907068C082C083907068E0F8A3E0FAC000C002907058E0F8A3E0FA84  
:201B3300C000C002D003D001D002D00012350ED083D082E8F0A3EAF0907068E0F8A3E0FABE  
:201B5300C000C00278007A10C000C002D003D001D002D00012352DE84A7003021B8590706B  
:201B730068C082C08378007A10D083D082E8F0A3EAF0907068C082C083907068E0F8A3E0B9  
:201B9300FAC000C00290706AE0F8A3E0FAC000C002D003D001D002D00012348ED083D082B6  
:201BB300E8F0A3EAF0907062C082C083907068E0F8A3E0FAC000C00290705EE0F8A3E0FAE4

00E090" 5EEFE950

:201BD300C000C002D003D001D002D000123512D083D082E8F0A3EAF090706AE0F8A3E0FA18  
:201BF300C000C00278017A00C000C002D003D001D002D00012342FE84A7003021C679070F6  
:201C130062E0F8A3E0FAC000C002907060E0F8A3E0FAC000C002D003D001D002D0001235B4  
:201C33002DE84A7003021C64907062C082C083907060E0F8A3E0FAD083D082E8F0A3EAF0A7  
:201C5300907085C082C08378C87A00D083D082E8F0021CBA907062E0F8A3E0FAC000C0021F  
:201C7300907060E0F8A3E0FAC000C002D003D001D002D000123553E84A7003021CBA9070BD  
:201C930062C082C083907060E0F8A3E0FAD083D082E8F0A3EAF0907085C082C08378C87AD7  
:201CB30000D083D082E8F0907016C082C083907084E0F87A00C000C002C000C002781E7A0F  
:201CD30000C000C002D003D001D002D000123508E84A70030230C278027A00C000C002D0FB  
:201CF30003D001D002D0001234CFC000C002D003D001D002D00012350EC000C00290705C4B  
:201D1300E0F8A3E0FAC000C00290705AE0F8A3E0FAC000C002D003D001D002D000123537E4  
:201D3300D083D082E8F0A3EAF0907070C082C08378007A00D083D082E8F0A3EAF09070763F  
:201D5300C082C08378007A00D083D082E8F0A3EAF090700AC082C08378137A00D083D082C6  
:201D7300E8F0A3EAF0907008C082C08378017A00D083D082E8F0A3EAF0907070C082C0838C  
:201D9300907070E0F8A3E0FAC000C002907016C082C083907076E0F8A3E0FAC000C002C041  
:201DB30000C002781E7A00C000C002D003D001D002D000123508E84A70030230C278027A9A  
:201DD30000C000C002D003D001D002D0001234CFC000C002D003D001D002D00012350EC036  
:201DF30000C002D003D00189828B831235A5C000C002D003D001D002D000123512D083D021  
:201E130082E8F0A3EAF0907076C082C083907076E0F8A3E0FAC000C002907008E0F8A3E02D  
:201E3300FAC000C002D003D001D002D000123512D083D082E8F0A3EAF0907076E0F8A3E0A9  
:201E5300FAC000C00290700AE0F8A3E0FAC000C002D003D001D002D0009070081235ABE8EA  
:201E73004A6003021D8C907070C082C083907070E0F8A3E0FAC000C002907072E0F8A3E0EE  
:201E9300FAC000C002D003D001D002D00012348EC000C00278027A00C000C002D003D001FD  
:201EB300D002D0001235B7D083D082E8F0A3EAF0907080E0F8A3E0FAC000C00278007A002C  
:201ED300C000C002D003D001D002D000123596E84A7003021F18907070C082C08390707007  
:201EF300E0F8A3E0FAC000C00278027A00C000C002D003D001D002D0001235B7D083D08299  
:201F1300E8F0A3EAF0907054E0F87A00C000C00278017A00C000C002D003D001D002D00076  
:201F330012342FC000C002907085E0F87A00C000C00278007A00C000C002D003D001D00254  
:201F5300D00012342FC000C00212359CC000C00278017A00C000C002D003D001D002D00087  
:201F730012354DC000C002D003D001D002D00012356EC000C002D003D001D002D00012352E  
:201F930075E84A70030225B907054E0F87A00C000C00278017A00C000C002D003D001D084  
:201FB30002D00012342FE84A7003021FEB907010C082C08378147A00D083D082E8F0A3EA71  
:201FD300F090700EC082C08378F67AFFD083D082E8F0A3EAF0022013907010C082C0837848  
:201FF300287A00D083D082E8F0A3EAF090700EC082C08378D87AFFD083D082E8F0A3EAF0DC  
:2020130090706CC082C083907062E0F8A3E0FAC000C00290705CE0F8A3E0FAC000C002D080  
:2020330003D001D002D000123537D083D082E8F0A3EAF0907068C082C08390705CE0F8A3DB  
:20205300E0FAC000C002907070E0F8A3E0FAC000C002D003D001D002D000123512C000C0AB  
:2020730002907062E0F8A3E0FAC000C002D003D001D002D000123537C000C00278027A00D8  
:20209300C000C002D003D001D002D0001235B7D083D082E8F0A3EAF090706CC082C08390EC  
:2020B300706CE0F8A3E0FAC000C002907068E0F8A3E0FAC000C002D003D001D002D00012C3  
:2020D3003537D083D082E8F0A3EAF090706CE0F8A3E0FAC000C002907010E0F8A3E0FAC01F  
:2020F30000C002D003D001D002D00012352DE84A700302212090706CC082C083907010E088  
:20211300F8A3E0FAD083D082E8F0A3EAF090706CE0F8A3E0FAC000C00290700EE0F8A3E091  
:20213300FAC000C002D003D001D002D000123553E84A700302216290706CC082C083907015  
:202153000EE0F8A3E0FAD083D082E8F0A3EAF090706EC082C083907062E0F8A3E0FAC000A5  
:20217300C00278117A00C000C002D003D001D002D0001235B7C000C00278467A00C000C087  
:2021930002D003D001D002D000123537C000C00290706CE0F8A3E0FAC000C002D003D001FD  
:2021B300D002D000123512D083D082E8F0A3EAF0907083C082C08378007A00D083D082E890  
:2021D300F0907083C082C083907083E0F87A00C000C00278017A00C000C002D003D001D0B4  
:2021F30002D0001235127A00D083D082E8F0907083E0F87A00C000C00290706EE0F8A3E08A  
:20221300FAC000C002D003D001D002D000123568E84A700302223012014D02223312015027  
:20223300907083E0F87A00C000C00278C87A00C000C002D003D001D002D000123582E84AB7  
:2022530070030221D40222A990705CE0F8A3E0FAC000C002907070E0F8A3E0FAC000C002BA  
:20227300D003D001D002D000123512C000C002907062E0F8A3E0FAC000C002D003D001D07D  
:2022930002D000123568E84A70030222A612014D0222A91201501225FF90705AC082C08396  
:2022B30090705CE0F8A3E0FAD083D082E8F0A3EAF0907084C082C083907084E0F87A00C0C1  
:2022D30000C00278017A00C000C002D003D001D002D0001235127A00D083D082E8F090701E

00E080-6E7E960

:2022F30084E0F87A00C000C00278137A00C000C002D003D001D002D00012352DE84A70038D  
:20231300022327907084C082C08378007A00D083D082E8F0221200FA1200FA22907087E023  
:20233300F87A00C000C00278017A00C000C002D003D001D002D00012342FE84A700302259A  
:20235300F1907057E0F87A00C000C00278017A00C000C002D003D001D002D00012342FE836  
:202373004A700302239090707CC082C08378017A00D083D082E8F0A3EAF00225F19070577B  
:20239300E0F87A00C000C00278027A00C000C002D003D001D002D00012342FE84A7003027E  
:2023B30023B81225F4907057E0F87A00C000C00278037A00C000C002D003D001D002D0001C  
:2023D30012342FE84A70030223F1907056C082C08378007A00D083D082E8F0907064C082C8  
:2023F300C083907089C082C083907074E0F8A3E0FAC000C00278087A00C000C002D003D00F  
:2024130001D002D00012348EC000C00278037A00C000C002D003D001D002D000123512C0DA  
:2024330000C002C000C00278E87A03C000C002D003D001D002D000123508E84A700302307A  
:20245300C278017A00C000C002D003D001D002D0001234CFC000C002D003D001D002D0000F  
:2024730012350EC000C002D003D00189828B83123528C000C00278207A00C000C002D0035D  
:20249300D001D002D00012348EC000C002780A7A00C000C002D003D001D002D0001235B79E  
:2024B300D083D082E8F0A3EAF0907064C082C083907064E0F8A3E0FAC000C00290706AE0A1  
:2024D300F8A3E0FAC000C002D003D001D002D00012348ED083D082E8F0A3EAF0907062C0BC  
:2024F30082C083907062E0F8A3E0FAC000C002907064E0F8A3E0FAC000C002D003D001D01C  
:2025130002D000123512D083D082E8F0A3EAF090706AE0F8A3E0FAC000C00278017A00C08F  
:2025330000C002D003D001D002D00012342FE84A700302259E907062E0F8A3E0FAC000C06A  
:2025530002907060E0F8A3E0FAC000C002D003D001D002D00012352DE84A700302259B907E  
:202573007062C082C083907060E0F8A3E0FAD083D082E8F0A3EAF0907057C082C837801ED  
:202593007A00D083D082E8F00225F1907062E0F8A3E0FAC000C002907060E0F8A3E0FAC06B  
:2025B30000C002D003D001D002D000123553E84A70030225F1907062C082C083907060E082  
:2025D300F8A3E0FAD083D082E8F0A3EAF0907057C082C08378017A00D083D082E8F00235F6  
:2025F300EFC29775F005D5F0FDD29722907086C082C08312361A750200D083D082E8F090D8  
:202613007086E0F87A00C000C00278087A00C000C002D003D001D002D00012342FE84A7004  
:2026330003022D1E907068C082C08378007A00D083D082E8F0A3EAF090707AC082C08378E7  
:20265300007A00D083D082E8F0A3EAF0907006C082C08378077A00D083D082E8F0A3EAF075  
:20267300907004C082C08378017A00D083D082E8F0A3EAF0907086C082C08312361A7502ED  
:2026930000D083D082E8F0907086E0F87A00C000C00278007A00C000C002D003D001D00266  
:2026B300D00012366CE84A7003022687907068C082C083907068E0F8A3E0FAC000C0029073  
:2026D3007086E0F87A00C000C002D003D001D002D000123512D083D082E8F0A3EAF0907084  
:2026F30089C082C08390707AE0F8A3E0FAC000C00278E87A03C000C002D003D00E4  
:2027130001D002D000123508E84A70030230C278017A00C000C002D003D001D002D000124E  
:2027330034CFC000C002D003D001D002D00012350EC000C002907086E0F87A00C000C0028A  
:2027530078207A00C000C002D003D001D002D0001235377A00D083D082E8F090707AC0825B  
:20277300C08390707AE0F8A3E0FAC000C002907004E0F8A3E0FAC000C002D003D001D00261  
:20279300D000123512D083D082E8F0A3EAF090707AE0F8A3E0FAC000C002907006E0F8A331  
:2027B300E0FAC000C002D003D001D002D0009070041235ABE84A6003022687907066C08282  
:2027D300C083907068E0F8A3E0FAC000C00278D27A00C000C002D003D001D002D000123591  
:2027F300B7D083D082E8F0A3EAF0907066C082C083907066E0F8A3E0FAC000C00278D27A29  
:2028130000C000C002D003D001D002D00012348ED083D082E8F0A3EAF0907068C082C08322  
:20283300907068E0F8A3E0FAC000C002907066E0F8A3E0FAC000C002D003D001D002D000C3  
:20285300123537D083D082E8F0A3EAF0907086C082C08312361A750200D083D082E8F090FC  
:202873007086E0F87A00C000C00278007A00C000C002D003D001D002D00012366CE84A706B  
:202893000302285F907086C082C083907086E0F87A00C000C00278207A00C000C002D003CD  
:2028B300D001D002D0001235377A00D083D082E8F0907086E0F87A00C000C002907068E07B  
:2028D300F8A3E0FAC000C002D003D001D002D00012366CE84A70030228FBD2179037B712B1  
:2028F30035881233D8123416907068E0F8A3E0FAC000C002907086E0F87A00C000C002D026  
:2029130003D001D002D00012342FE84A7003022637D2179037B91235881233D81234169074  
:2029330070778C082C08378017A00D083D082E8F0A3EAF0907068C082C08378007A00D083F8  
:20295300D082E8F0A3EAF090707AC082C08378007A00D083D082E8F0A3EAF0907002C0828E  
:20297300C08378077A00D083D082E8F0A3EAF0907000C082C08378017A00D083D082E8F019  
:20299300A3EAF0907086C082C08312361A750200D083D082E8F0907086E0F87A00C000C0EE  
:2029B3000278007A00C000C002D00





20313300E40529B52902052AF52830113820133520120EE52BB52A07E52CB5290280025064  
203153002420171F20081CE5C620E217D0E0D213C212C211C082C083903173C082C083322E  
203173000230FDD212D0E0D0D032E4F529F52AF528F52CF52BC211C212C213228528E07588  
20319300F005A485E00085F002227501057503001231AA85002822C3E94B7005D30232065D  
2031B300227F007E007C007D0075F0101231EFEE33FEEF33FFC3EF9B40117005C3EE994060  
2031D3000AC3EE99FEEF9BFFD38001C31231F7D5F0DA8C008D028F038E01C322C3E833F81A  
2031F300EA33FA22EC33FCED33FD229037F374140237399037E4740A023739903810742866  
20321300023739589F0D21043A882D28C22907489E8F0A3EAF0A4A88222E859F822E84991  
20323300F82290748BE8F0A3EAF0D219D2BCC2BAD2B943A884229074F87475F0A37400F090  
203253009074F67475F0A37400F0751000751101C208C21722C2ACC082C083F89074F6E0EB  
20327300F9A3E0F5828983E8F0A3E582B45003907500A882E5829074F7F088828983A3E534  
2032930082B45003907500A8829074F8A3E0B50008D2AC0000C2AC80F4200804D208D29956  
2032B300D083D082D2AC22209846C299300E03023390C0E0C082C083C000C001C002C0032C  
2032D3009074F7E0FA9074F8E0F9A3E0FBB50205C208023382850382850183E0F599A3E56D  
2032F30082B45003907500E5829074F9F0023382C0E0E599C29820150302338EB4030620D2  
203313000F030233EEB41304D20E806FB41109300E69D299C20E8063C082C083C000C00132  
203333300C002C003FA90749FE0F890749DAB829074A0E0F9A3E0F5828983A3E582B503036A  
2033530090748DEAF0E583A9829074A0F0E9A3F0B5001C90749EE0F5838882A3E582B503B5  
203373000390748DE583A98290749FE0E9A3F0D003D002D001D000D083D082D0E0D0D03298  
20339300E060061233A4A380F72008F222008FD221233AB123268227911B42000400207DE  
2033B30022B40D03770122B40081722B40A07D0E0D0E00233D22C000C0E0123268D0E06E  
2033D300D000D8F322740D1233A470A123268E5106005F8E41233CA22C2AF9033F6C082B6  
2033F300C08332C3721B721A5008903405C082C083323008033099FD851581301D04C2AFB3  
2034130080FE222008FDC21712341F22301A0C90342EC082C083C0D00230D722123435026F  
203433003462C2D5EB33EA123459903441738012D322200EACBFAE8C9F812347620D5017A  
20345300B32212347622333333C3540622B2D5A2D58008B3B0D58003B372D5400302348595  
20347300023489C2D5EAB50308C2D5E8B50102D2D522E4F8FA2278FF7AFF22EB33EA1234E1  
2034930059C0E090349A73800F802780031234F61234FFEA5B33401F1234CFEA30E70A33EB  
2034B300487014D0E030D01122D0E030D0031234F6221234F680E115810231FEC3E889F0B1  
2034D300A4FCADF0E960088AF0A420D2122DFDEB600688F0A420D2072D4004FAA80422D38D  
2034F3000231FEE4C398F8E49AFAC322E4C399F9E49BFBC322123476023467D2D58002C218  
20351300D5C3E829F8EA3BFA10D504A2D2C250030231FE22E0F87A0022123435023466B6  
20353300D2D58002C2D5C3E899F8EA9BFA10D504A2D2C250030231FE221234760234660E  
20355300123435023467E4A30E70A1234F6EA30E70302320E2212343502346CE859F8EA5BF6  
20357300FA22E849F8EA4BFA2212347C02346612343502346BE49360061233A4A380F62029  
2035930008FD22123476023462E52AA829B52AF9FA22E0F8A3E0FA22C3A3E033400302355F  
2035B30068023582EB33EA1234599035C17380118013801E1234F61234FF1231AA1234FFC2  
2035D300221231AA221234F61231AA1234F61234FF221234FF1231AA1234F622D007D0069E  
2035F300D005D004D003D002D001D000D0F0D0E0D0D020D509D083D082D0D0C21B32D08314  
20361300D082D0D0C21322C0A853A8EB9074A1E0F890749DAB8290749EE0F9A3E0B500045E  
203633007800801A8983F582A3E582B5030390748DE0F8AB82E58390749EF0EBA3F0D0A897  
203653002290748DAA83A98290749EEAF0A3E9F09074A0EAF0A3E9F022123476023460C224  
20367300111233A080FE7433C3301B02940EF5F0E581B51500400C7477F42581400525F022  
20369300400122851581023269C903871746402739C082C0831233D890389712358812362A  
2036B300BFD083D0821235881233D82290389E12358890755B1235A5907551A982AB831243  
2036D30036E29075511233939038A512358822C3E84A7006C0E07C01800E7C00790A1237D5  
2036F30007C0010CE84A70F47930D0E029F0A3DCF9E4F022AF02AE00AB0178077902EF8BEE  
20371300F084AAF0FF18E6D7C48BF084AAF0C4C6C4D7C48BF084AAF0D6A806AA07A9F022E5  
20373300E8F0A3EAF022B42800500122C0E0C000C002C0E0C082C0839075591235A590751A  
203753005B123733D083D082D0E0B4280040061236A42008FDD002D000D0E090755DF0B49F  
20377300500122B42800500122B4960122301B0B301A0890378EC082C08332C2AF80FEC2A2  
20379300AF301B0B301A089037A2C082C08332C2AF80FE5045533A54454D5020322E34356F  
2037B300610061004200470047004200417272617920496E664657820546F6F204269670087  
2037D300537472696E67205472756E636174656400446976696465206279205A65726F0085  
2037F3004D617468204F766572666C6F77004D6174682055566E646572666C6F77004D617436  
203813006820417267756D665E742042616400496E7075742056616C756520546F6F2042C2  
20383



```
:20385300204E656761746976650049737461636B2F42353220436F6E666C696374004973BD
:203873007461636B204572726F7200566563746F722072656C6F636174696F6E20464149B5
:203893004C4544004552524F522000284C696E652000293A2000535441525455502F52454A
:0438B3005345540025
:00003301CC
```

Stepper processor(IC 3) hex file(Step23.hex):

```
:03000000020033C8
:05000B00C0D0021EB68A
:05001300C0D0021E68D0
:05002300C0D0021FDF48
:05002B00C0D0021C3CE6
:20003300C2AFC2D4C2D3752100752200752300D215D28AAF81758130858115C0D0C0F0C098
:20005300E0C08753C6FB121C837527405388CF5389F0121F0E78FC121C57121F6E1221AC99
:20007300C20EC20FC208759850D2ACD2AFE587D0E0D000D001D0D030E20530E4028012206A
:20009300E602800DBF0702800890245374FA0222E475C7AA75C755438740907058C082C030
:2000B3008378177A00C000C002902351C082C083D003D001D005D004D002D000121C8F905A
:2000D3007024C082C08378197A00C000C002C000C002781A7A00C000C002D003D001D002E1
:2000F300D000121CCCE84A7003021CBC78027A00C000C002D003D001D002D000121D00C0F9
:2001130000C002D003D001D002D000121D3FC000C00278547A3DC000C002D002D000D083DA
:20013300D082E8F0A3EAF0907009C082C083783C7A00D083D082E8F002068D907004C082F1
:20015300C08378007A00D083D082E8F0907010E0F8A3E0FAC000C00278807A3EC000C002C1
:20017300D003D001D002D000121D59E84A70030201B4907010C082C083907010E0F8A3E042
:20019300FAC000C00278807A3EC000C002D003D001D002D000121D8ED083D082E8F0A3EA91
:2001B300F0907010E0F8A3E0FAC000C00278007A00C000C002D003D001D002D000121DA4C8
:2001D300E84A7003020209907010C082C083907010E0F8A3E0FAC000C00278807A3EC0006E
:2001F300C002D003D001D002D000121D43D083D082E8F0A3EAF0907010E0F8A3E0FAC00053
:20021300C00278407A1FC000C002D003D001D002D000121D59E84A7003020272907004C089
:2002330082C08378017A00D083D082E8F0120631907010C082C08378807A3EC000C00290D6
:200253007010E0F8A3E0FAC000C002D003D001D002D000121D8ED083D082E8F0A3EAF090A7
:200273007012C082C083907010E0F8A3E0FAC000C00290700EE0F8A3E0FAC000C002D003C5
:20029300D001D002D000121DAAD083D082E8F0A3EAF090700CC082C08390700CE0F8A3E00D
:2002B300FAC000C002907012E0F8A3E0FAC000C002D003D001D002D000121D43D083D08269
:2002D300E8F0A3EAF090700CE0F8A3E0FAC000C00278807A3EC000C002D003D001D002D05B
:2002F30000121DFDE84A700302032D90700CC082C08390700CE0F8A3E0FAC000C00278807C
:200313007A3EC000C002D003D001D002D000121D8ED083D082E8F0A3EAF0907010E0F8A308
:20033300E0FAC000C00278017A00C000C002D003D001D002D000121D59E84A70030205D28D
:20035300907015E0F87A00C000C00278017A00C000C002D003D001D002D000121DA4E84AE1
:200373007003020389907015C082C08378017A00D083D082E8F0907015E0F87A00C000C078
:2003930002780B7A00C000C002D003D001D002D000121D59E84A70030203BF907015C0823B
:2003B300C083780B7A00D083D082E8F090701EC082C083907010E0F8A3E0FAC000C002786B
:2003D300027A00C000C002D003D001D002D000121E03D083D082E8F0A3EAF090701EE0F8A3
:2003F300A3E0FAC000C002907015E0F87A00C000C002D003D001D002D000121E57E84A7093
:200413000302042E90701EC082C083907015E0F87A00D083D082E8F0A3EAF090701CC08230
:20043300C08378017A00D083D082E8F0A3EAF090701AC082C08378157A00D083D082E8F056
:20045300A3EAF0907012C082C083907010E0F8A3E0FAC000C00290701EE0F8A3E0FAC0005B
:20047300C00290701CE0F8A3E0FAC000C002D003D001D002D000121D00C000C00278027AC9
:2004930000C000C002D003D001D002D000121D00C000C002D003D001D002D000121D8AD001
:2004B30083D082E8F0A3EAF0907020C082C083907012E0F8A3E0FAC000C00278C87A00C0F7
:2004D30000C002D003D001D002D000121E03D083D082E8F0A3EAF0907022C082C08390702D
:2004F30012E0F8A3E0FAC000C00278C87A00C000C002907020E0F8A3E0FAC000C002D003FA
:20051300D001D002D000121D00C000C002D003D001D002D000121D8AD083D082E8F0A3EA9B
:20053300F0907022E0FE907020E0FF90701EE0F890701CE0F990701AE0FA74FF88F0C0F03F
:200553009A89F0C0F0D2A0F5F0D5F0FDC2A0F5F0D5F0FDD0F0D5F0EBD0F0D5F0E1FBEF24BF
:2005730001FFDF0302059475F0C8C0F0D2A08BF0D5F0FDC2A08BF0D5F0FDD0F0D5F0EBDF71
```

20059300E6EE2401FEDE032025ADD2A08BF0D5F0FDC2A08BF0D5F0FDDFE0EB88F0C0F02A63  
2005B30089F0C0F0D2A0F5F0D5F0FDC2A0F5F0D5F0FDD0F0D5F0EBD0F0D5F0E10206089067  
2005D3007010E0F8A3E0FAC000C00278017A00C000C002D003D001D002D000121E5DE84A37  
2005F3007003020608D2A075F0FFD5F0FDC2A075F0FFD5F0FDD907004E0F87A00C000C0026D  
2006130078017A00C000C002D003D001D002D000121E5DE84A700302063012063122907037  
200633000EE0F8A3E0FAC000C00278FF7AFFC000C002D003D001D002D000121E5DE84A70DB  
200653000302067090700EC082C08378017A00D083D082E8F0A3EAF0D2A202068690700E7C  
20067300C082C08378FF7AFFD083D082E8F0A3EAF0C2A222D2A122C2A12278007A00121E36  
200693006378007A00121F2E121F3BA8887A00C000C00278F77A00C000C002D003D001D01C  
2006B30002D000121F528888A8887A00C000C00278407A00C000C002D003D001D002D0009C  
2006D300121F5688889019BAA882AA83121F5AD297C2A3D2A2C4A1D2179070581220AD  
2006F300B81220FD1213B90700EC082C08378017A00C000C002D002D000D083D082E8F06B  
20071300A3EAF0907016C082C08378007A00D083D082E8F0907008C082C08378007A00D0F0  
2007330083D082E8F0907018C082C08378007A00D083D082E8F0907019C082C08378197A44  
2007530000D083D082E8F078007A00121E6378007A00121F2E121F20C000C00278FA7A0074  
20077300C000C002D003D001D002D000121E57E84A700302078C121AFE120C9E907016E001  
20079300F87A00C000C00278017A00C000C002D003D001D002D000121E5DE84A700302075C  
2007B30068907018E0F87A00C000C00278017A00C000C002D003D001D002D000121D59E8A7  
2007D3004A7003020B2590700AC082C08378017A00D083D082E8F0A3EAF0907002C082C097  
2007F30083907018E0F87A00D083D082E8F0A3EAF0907000C082C08378017A00D083D082E2  
20081300E8F0A3EAF090700AE0F8A3E0FAC000C002907019E0F87A00C000C002D003D001FE  
20083300D002D000122154E84A700302096E907010C082C083907024C082C08390700AE036  
20085300F8A3E0FAC000C002C000C002781A7A00C000C002D003D001D002D000121CCCE856  
200873004A7003021CBC78027A00C000C002D003D001D002D000121D00C000C002D003D0BE  
2008930001D002D000121D3FC000C002D003D00189828B8312215AC000C002907024C08280  
2008B300C083907019E0F87A00C000C002C000C002781A7A00C000C002D003D001D002D09F  
2008D30000121CCCE84A7003021CBC78027A00C000C002D003D001D002D000121D00C000E1  
2008F300C002D003D001D002D000121D3FC000C002D003D00189828B8312215AC000C00221  
20091300D003D001D002D000121D8AD083D082E8F0A3EAF0907010C082C083907010E0F84E  
20093300A3E0FAC000C00290700EE0F8A3E0FAC000C002D003D001D002D000121DAAD0834E  
20095300D082E8F0A3EAF0907015C082C083780A7A00D083D082E8F012014E907019C0820E  
20097300C08390700AE0F8A3E0FA7A00D083D082E8F0121A2090700AC082C08390700AE006  
20099300F8A3E0FAC000C002907000E0F8A3E0FAC000C002D003D001D002D000121D43D0EE  
2009B30083D082E8F0A3EAF090700AE0F8A3E0FAC000C002907002E0F8A3E0FAC000C00240  
2009D300D003D001D002D000907000122160E84A6003020818907010C082C083907024C0FB  
2009F30082C08378017A00C000C002C000C002781A7A00C000C002D003D001D002D0001242  
200A13001CCCE84A7003021CBC78027A00C000C002D003D001D002D000121D00C000C002EF  
200A3300D003D001D002D000121D3FC000C002D003D00189828B8312215AC000C0029070A1  
200A530024C082C083907019E0F87A00C000C002C000C002781A7A00C000C002D003D00139  
200A7300D002D000121CCCE84A7003021CBC78027A00C000C002D003D001D002D000121D5D  
200A930000C000C002D003D001D002D000121D3FC000C002D003D00189828B8312215AC081  
200AB30000C002D003D001D002D000121D8AD083D082E8F0A3EAF0907010C082C0839070D3  
200AD30010E0F8A3E0FAC000C00290700EE0F8A3E0FAC000C002D003D001D002D000121DC2  
200AF300AAD083D082E8F0A3EAF0907015C082C083780A7A00D083D082E8F012014E9070CB  
200B130019C082C08378017A00D083D082E8F002C87907019E0F87A00C000C00278017A3F  
200B330000C000C002D003D001D002D000122154E84A7003020C73907010C082C083907098  
200B530024C082C08378017A00C000C002C000C002781A7A00C000C002D003D001D002D00E  
200B730000121CCCE84A7003021CBC78027A00C000C002D003D001D002D000121D00C0003E  
200B9300C002D003D001D002D000121D3FC000C002D003D00189828B8312215AC000C0027E  
200BB300907024C082C083907019E0F87A00C000C002C000C002781A7A00C000C002D003A9  
200BD300D001D002D000121CCCE84A7003021CBC78027A00C000C002D003D001D002D0005A  
200BF300121D00C000C002D003D001D002D000121D3FC000C002D003D00189828B8312210B  
200C13005AC000C002D003D001D002D000121D8AD083D082E8F0A3EAF0907010C082C08357  
200C3300907010E0F8A3E0FAC000C00290700EE0F8A3E0FAC000C002D003D001D002D0008F  
200C5300121DAAD083D082E8F0A3EAF0907015C082C083780A7A00D083D082E8F012014E3A  
200C7300907019C082C08378017A00D083D082E8F0121A20907016C0

:200CB30017E0F87A00C000C00278097A00C000C002D003D001D002D000121E5DE84A700341  
 :200CD300021834120687907008C082C08378007A00D083D082E8F090700EE0F8A3E0FAC0F5  
 :200CF30000C00278FF7AFFC000C002D003D001D002D000121E5DE84A7003020D13120631CA  
 :200D130090700CC082C08378007A00D083D082E8F0A3EAF0907015C082C08378057A00D0E2  
 :200D330083D082E8F0907005C082C08378007A00D083D082E8F0907017C082C083122172B9  
 :200D5300750200D083D082E8F0907017E0F87A00C000C00278007A00C000C002D003D00189  
 :200D7300D002D000122154E84A7003020D49907017E0F87A00C000C00278527A00C000C08B  
 :200D930002D003D001D002D000121E5DE84A7003020DA9120631907017E0F87A00C000C0DC  
 :200DB3000278507A00C000C002D003D001D002D000121E5DE84A7003020DE5907010C0829C  
 :200DD300C08378017A00D083D082E8F0A3EAF012014E907017E0F87A00C000C00278447A4E  
 :200DF30000C000C002D003D001D002D000121E5DE84A7003020E58907010C082C083780A67  
 :200E13007A00D083D082E8F0A3EAF012014E90700CE0F8A3E0FAC000C002786B7A3EC000AC  
 :200E3300C002D003D001D002D000121D59E84A7003020E58907005C082C08378017A00D0B5  
 :200E530083D082E8F0907017E0F87A00C000C00278467A00C000C002D003D001D002D000E7  
 :200E7300121E5DE84A7003020ECB907010C082C08378147A00D083D082E8F0A3EAF01201AA  
 :200E93004E90700CE0F8A3E0FAC000C002786B7A3EC000C002D003D001D002D000121D5923  
 :200EB300E84A7003020ECB907005C082C08378017A00D083D082E8F0907017E0F87A00C07C  
 :200ED30000C00278477A00C000C002D003D001D002D000121E5DC000C002907005E0F87AD6  
 :200EF30000C000C00278007A00C000C002D003D001D002D000122154C000C002D003D001F6  
 :200F1300D002D0001221D2E84A7003020D49907005E0F87A00C000C00278007A00C000C0CF  
 :200F330002D003D001D002D000121E5DE84A700302181590700EE0F8A3E0FAC000C0027898  
 :200F5300017A00C000C002D003D001D002D000121E5DE84A7003020F6F120631907017C069  
 :200F730082C083122172750200D083D082E8F0907017E0F87A00C000C00278007A00C00063  
 :200F9300C002D003D001D002D000122154E84A7003020F6F907017C082C083907017E0F8FF  
 :200FB3007A00C000C00278207A00C000C002D003D001D002D000121D8E7A00D083D082E884  
 :200FD300F0907010C082C08378967A00C000C002907017E0F87A00C000C002D003D001D010  
 :200FF30002D000121D8ED083D082E8F0A3EAF0907015C082C083780A7A00D083D082E8F042  
 :2010130012014E907024C082C08378017A00C000C002C000C002781A7A00C000C002D0035B  
 :20103300D001D002D000121CCCE84A7003021CBC78027A00C000C002D003D001D002D000F5  
 :20105300121D00C000C002D003D001D002D000121D3FC000C00278007A00D083D082E8F027  
 :20107300A3EAF0907007C082C08378017A00D083D082E8F0907010C082C08378367A05D052  
 :2010930083D082E8F0A3EAF0907015C082C083780A7A00D083D082E8F012014E90700CC0D3  
 :2010B30082C08378CA7AFAD083D082E8F0A3EAF0120631907017C082C083122172750200A7  
 :2010D300D083D082E8F0907017E0F87A00C000C00278007A00C000C002D003D001D002D0DB  
 :2010F30000122154E84A70030210C6907017E0F87A00C000C002784E7A00C000C002D00359  
 :20111300D001D002D000121E5DE84A7003021178907007C082C083907007E0F87A00C000E7  
 :20113300C00278017A00C000C002D003D001D002D000121D437A00D083D082E8F0907010A6  
 :20115300C082C083789B7A02D083D082E8F0A3EAF0907015C082C083780A7A00D083D08233  
 :20117300E8F012014E907017E0F87A00C000C00278477A00C000C002D003D001D002D00037  
 :20119300121E5DE84A70030210C6907017C082C083122172750200D083D082E8F0907017E6  
 :2011B300E0F87A00C000C00278007A00C000C002D003D001D002D000122154E84A70030260  
 :2011D300119D907007C082C083907017E0F87A00C000C00278207A00C000C002D003D0019F  
 :2011F300D002D000121D8E7A00D083D082E8F0907010C082C08378807A3EC000C0029070BF  
 :20121300CE0F8A3E0FAC000C002D003D001D002D000121D8ED083D082E8F0A3EAF09070DB  
 :2012330015C082C083780A7A00D083D082E8F012014E90700EE0F8A3E0FAC000C00278FFCB  
 :201253007AFFC000C002D003D001D002D000121E5DE84A700302126E120631D21790235F42  
 :201273001221C41220FD12213B907017C082C083122172750200D083D082E8F0907017E09B  
 :20129300F87A00C000C00278007A00C000C002D003D001D002D000122154E84A700302124D  
 :2012B3007C907007E0F87A00C000C00278017A00C000C002D003D001D002D000121E5DE894  
 :2012D3004A700302130C907018C082C08378017A00D083D082E8F0D2179023611221C4120A  
 :2012F30020FD12213B907019C082C08378017A00D083D082E8F0021812907014C082C0837D  
 :2013130078007A00D083D082E8F0907018C082C08378017A00D083D082E8F0907017C082E5  
 :20133300C083122172750200D083D082E8F0907017E0F87A00C000C00278007A00C000C061  
 :2013530002D003D001D002D000122154E84A700302132E907017E0F87A00C000C00278520E  
 :201373007A00C000C002D003D001D002D000121E5DE84A700302138E120631907017E0F80B  
 :201393007A00C000C00278507A00C000C002D003D001D002D000121E5DE84A70030213CA23  
 :2013B300907010C082C08378017A00D083D082E8F0A3EAF012014E907017E0F87A00C0000E

:2013D300C00278447A00C000C002D003D001D002D000121E5DE84A7003021417907010C00B  
 :2013F30082C083780A7A00D083D082E8F0A3EAF0907015C082C08378047A00D083D082E802  
 :20141300F012014E907017E0F87A00C000C00278467A00C000C002D003D001D002D000126B  
 :201433001E5DE84A70030214E907010C082C08378147A00D083D082E8F0A3EAF0907015D9  
 :20145300C082C08378087A00D083D082E8F090700CE0F8A3E0FAC000C00278287A3CC00084  
 :20147300C002D003D001D002D000121DA4C000C002907018E0F87A00C000C002907007E029  
 :20149300F87A00C000C002D003D001D002D0001221D9C000C002D003D001D002D0001221F8  
 :2014B300DFE84A70030214C112014E0214E0D2179023631221C41220FD12213B907014C000  
 :2014D30082C08378017A00D083D082E8F0907017E0F87A00C000C00278547A00C000C00211  
 :2014F300D003D001D002D000121E5DE84A700302152D907010C082C083782C7A01D083D046  
 :2015130082E8F0A3EAF0907015C082C083780A7A00D083D082E8F012014E907017E0F87A04  
 :201533000C000C00278477A00C000C002D003D001D002D000121E5DE84A70030216C39078  
 :201553007017C082C083122172750200D083D082E8F0907017E0F87A00C000C00278007AF6  
 :201573000C000C002D003D001D002D000122154E84A7003021552907017C082C08390705F  
 :2015930017E0F87A00C000C00278207A00C000C002D003D001D002D000121D8E7A00D083E9  
 :2015B300D082E8F0907018C082C083907018E0F87A00C000C00278017A00C000C002D0031D  
 :2015D300D001D002D000121D437A00D083D082E8F0907024C082C083907018E0F87A00C049  
 :2015F3000C002C000C002781A7A00C000C002D003D001D002D000121CCCE84A7003021C03  
 :20161300BC78027A00C000C002D003D001D002D000121D00C000C002D003D001D002D00048  
 :20163300121D3FC000C00290700CE0F8A3E0FAC000C00278957A00C000C002D003D001D047  
 :2016530002D000121D8AC000C002907017E0F87A00C000C002D003D001D002D000121D3FCB  
 :20167300D083D082E8F0A3EAF0907018E0F87A00C000C002907007E0F87A00C000C002D0C6  
 :2016930003D001D002D0001221E6E84A70030216C3D2179023651221C41220FD12213B9003  
 :2016B3007014C082C08378017A00D083D082E8F0907014E0F87A00C000C00278017A00C003  
 :2016D3000C002D003D001D002D000121E5DE84A700302132E907010C082C08378807A3E35  
 :2016F300C000C00290700CE0F8A3E0FAC000C002D003D001D002D000121D8ED083D082E8E2  
 :20171300F0A3EAF0907010E0F8A3E0FAC000C00278007A00C000C002D003D001D002D000A8  
 :20173300121DA4E84A700302176C907010C082C083907010E0F8A3E0FAC000C00278807AAB  
 :201753003EC000C002D003D001D002D000121D43D083D082E8F0A3EAF0907010E0F8A3E099  
 :20177300FAC000C00278407A1FC000C002D003D001D002D000121D59E84A70030217C412A5  
 :201793000631907010C082C08378807A3EC000C002907010E0F8A3E0FAC000C002D003D0AE  
 :2017B30001D002D000121D8ED083D082E8F0A3EAF0907015C082C083780A7A00D083D08281  
 :2017D300E8F012014E90700EE0F8A3E0FAC000C00278FF7AFFC000C002D003D001D002D020  
 :2017F30000121E5DE84A7003021801120631907019C082C08378017A00D083D082E8F00230  
 :201813001823D2179023671221C41220FD12213B907017C082C08378007A00D083D082E8C8  
 :20183300F0907017E0F87A00C000C00278087A00C000C002D003D001D002D000121E5DE883  
 :201853004A700302186A907017C082C08378007A00D083D082E8F0907017E0F87A00C00000  
 :20187300C002780A7A00C000C002D003D001D002D000121E5DE84A70030218A0907017C00C  
 :2018930082C08378007A00D083D082E8F0907017E0F87A00C000C002780C7A00C000C00296  
 :2018B300D003D001D002D000121E5DE84A70030218D6907017C082C08378007A00D083D0FC  
 :2018D30082E8F0907017E0F87A00C000C002780D7A00C000C002D003D001D002D000121EB9  
 :2018F3005DE84A7003021972907017C082C083122172750200D083D082E8F0907017E0F822  
 :201913007A00C000C00278007A00C000C002D003D001D002D000122154E84A70030218FBBBD  
 :20193300907017E0F87A00C000C00278217A00C000C002D003D001D002D000121E5DE84A0F  
 :20195300700302195E12068702196112068A907017C082C08378007A00D083D082E8F09030  
 :201973007017E0F87A00C000C002780E7A00C000C002D003D001D002D000121E5DE84A7002  
 :20199300030219A8907017C082C08378007A00D083D082E8F0907017C082C08378007A00D5  
 :2019B300D083D082E8F022907016E0F87A00C000C00278007A00C000C002D003D001D002A1  
 :2019D300D000121E5DE84A7003021A1D907018E0F87A00C000C00278007A00C000C002D089  
 :2019F30003D001D002D000121D59E84A7003021A1AD2A4907016C082C08378017A00D083A4  
 :201A1300D082E8F0021A1D121AF30221EC78007A00121F2E121F20C000C002907009E0F81D  
 :201A33007A00C000C002D003D001D002D000122217E84A7003021A27D2A378007A00121F86  
 :201A53002E907006C082C083A8A07A00C000C00278207A00C000C002D003D001D002D0009C  
 :201A7300121F52D083D082E8F0907006E0F87A00C000C00278007A00C000C002D003D00161  
 :201A9300D002D000121E5DE84A7003021A54C2A3C2A4907006C082C083A8A07A00C000C057  
 :201AB3000278207A00C000C002D003D001D002D000121F52D083D082E8F0907006E0F87ADF  
 :201AD30000C000C00278007A00C000C002D003D001D002D000122154E84A7003021AA522A8

:201AF300C29775F005D5F0FDD29722907010C082C083907024C082C08378197A00C000C0FA  
 :201B130002C000C002781A7A00C000C002D003D001D002D000121CCCE84A7003021CBC7869  
 :201B3300027A00C000C002D003D001D002D000121D00C000C002D003D001D002D000121D28  
 :201B53003FC000C002D003D00189828B8312215AC000C002907024C082C083907019E0F84B  
 :201B73007A00C000C002C000C002781A7A00C000C002D003D001D002D000121CCCE84A7064  
 :201B930003021CBC78027A00C000C002D003D001D002D000121D00C000C002D003D001D074  
 :201BB30002D000121D3FC000C002D003D00189828B8312215AC000C002D003D001D002D03E  
 :201BD30000121D8AD083D082E8F0A3EAF0907010C082C083907010E0F8A3E0FAC000C002C3  
 :201BF30090700E0F8A3E0FAC000C002D003D001D002D000121DAAD083D082E8F0A3EAF0D4  
 :201C1300907015CE082C083780A7A00D083D082E8F012014E907019C082C08378197A00D0C4  
 :201C330083D082E8F02202221D8587E020E5FAD0D0D0E0D0E0D2D574AA75F055900033C094  
 :201C530082C0833253883F53892F438920438780888D758B00D28E2275C7AA75C755D2BF80  
 :201C730022121C6B75C7AA75C755438704D21D2275C7AA75C7555387FBC21D2289828B83DA  
 :201C9300E493C082C08388828A83F06017DC09E4F015811581021CC4A3A882AA83D083D0D3  
 :201CB30082A380DC158115812290236974460222E490237D74500222E4121CE8021CD9B2A9  
 :201CD300D5A2D58008B3B0D58003B372D54003021CF7021CFBC2D5EAB50308C2D5E8B50181  
 :201CF30002D2D522E4F8FA2278FF7AFF22C3E889F0A4FCADF0E960088AF0A420D2122DFDFE  
 :201D1300EB600688F0A420D2072D4004FAA80422D3021D2790239D74140222E490238E7463  
 :201D33000A0222E49023BA74280222E4D2D58002C2D5C3E829F8EA3BFA10D504A2D2C2D2D1  
 :201D53005003021D2722121D5F021CD8C2D5EB33EA121D83901D6B738012D3222200EACBF7  
 :201D7300FAE8C9F8121CE820D501B322121CE82233333C35406222D2D5800C2D2D5C3E899B8  
 :201D9300F8E49BFA10D504A2D2C2D25003021D2722121D5F021CD9EB33EA121D83C0E0909E  
 :201DB3001DB67380F802780D121DEB121DF4EA5B33401F121D00EA30E70A33487014D09F4  
 :201DD300E030D01122D0E030D003121DEB22121DEB80E11581021D27E4C398F8E49AFAC325  
 :201DF30022E4C399F9E49BFBC322121D5F021CDDC3E94B7005D3021D2F227F007E007C0065  
 :201E13007D0075F010121E48EE33FEF33FFC3EF9B40117005C3EE99400AC3EE99FEF9B8C  
 :201E3300FFD38001C3121E50D5F0DA8C008D028F038E01C322C3E833F8EA33FA22EC33FC0F  
 :201E5300ED33FD22121CE8021CD8121D5F021CD488298A2A2230191C301705D21AD0D022BE  
 :201E73002008FDC082C083907071C2D5D21BC21AC0D0800ED0D032C082C08390706FD2D549  
 :201E9300C0D0C0E0C0F0C000C001C002C003C004C005C006C007D215E0F8A3E0F5838500F4  
 :201EB30082E473301052C0E085278CE52804B4C808E40529B52902052AF52830113820134D  
 :201ED3003520120EE52BB52A07E52CB529028002502420171F20081CE5C620E217D0E0D2BD  
 :201EF30013C212C211C082C083901F04C082C08332021E8ED212D00E0D032E4F529F52AF1  
 :201F1300F528F52CF5B2C211C212C213228528E075F005A485E00085F00222750105750326  
 :201F330000121E03850028225389F0D21043A882D28C2290706FE8F0A3EAF043A88222E8B6  
 :201F530059F822E849F822907071E8F0A3EAF0D219D2BCC2BAD2B943A884229070DE747116  
 :201F7300F0A37400F09070DC7471F0A37400F0751000751101C208C21722C2ACC082C083DB  
 :201F9300F89070DCE0F9A3E0F5828983E8F0A3E582B45003907100A882E5829070DDF088AB  
 :201FB300828983A3E582B45003907100A8829070DEA3E0B50008D2AC0000C2AC80F420089E  
 :201FD30004D208D299D083D082D2AC22209846C299300E030220B5C0E0C082C083C000C04A  
 :201FF30001C002C0039070DDE0FA9070DEE0F9A3E0FBB50205C2080220A7850382850183FA  
 :20201300E0F599A3E582B45003907100E5829070DFF00220A7C0E0E599C29820150302205C  
 :20203300B3B400306200F03022113B41304D20E806FB41109300E69D299C20E8063C082C086  
 :2020530083C000C001C002C003FA907085E0F8907083AB82907086E0F9A3E0F5828983A3D5  
 :20207300E582B50303907073EAF0E583A982907086F0E9A3F0B5001C907084E0F5838882A2  
 :20209300A3E582B50303907073E583A982907084F0E9A3F0D003D002D001D000D083D08287  
 :2020B300D0E0D0D032E060061220C9A380F72008FD222008FD221220D0121F8D227911B482  
 :2020D300200040020722B40D03770122B408021722B40A07D0E0D0E002210222C000C0E041  
 :2020F300121F8DD0E0D000D8F322740D1220C9740A121F8DE5106005F8E41220EF22C2AF00  
 :2021130090211BC082C08332C3721B721A500890212AC082C083323008033099FD85158147  
 :20213300301D04C2AF80FE222008FDC21712214422301AOC902153C082C083C0D0021E689C  
 :2021530022121CE8021CD2E0F8A3E0FA22C3A3E033400302216C021DFD121D5F021CDEC01C  
 :20217300A853A8EB907087E0F8907083AB82907084E0F9A3E0B500047800801A8983F58281  
 :20219300A3E582B50303907073E0F8AB82E583907084F0E9A3F0D0A822907073AA83A98230  
 :2021B3009

```
:02221300D0C21322121CE8021CDDC2111220C580FE7433C3301B02940EF5F0E581B515001D
:02223300400C747FF42581400525F040012285158102224790241B74640222E4C082C0833B
:022253001220FD9024411221C412226AD083D0821221C41220FD229024481221C490715B76
:0222730012215A907151A982AB8312228D9071511220B890244F1221C422C3E84A7006C0CF
:02229300E07C01800E7C00790A1222B2C0010CE84A70F47930D0E029F0A3DCF9E4F022AF69
:0222B30002AE00AB0178077902EF8BF084AAF0FF18E6D7C48BF084AAF0C4C6C4D7C48BF098
:0222D30084AAF0D6A806AA07A90702E8F0A3EAF022B42800500122C0E0C000C002C0E0C095
:0222F30082C08390715912215A90715B1222DED083D082D0E0B428004000612224F2008FD92
:0231300D002D000D0E090715DF0B4500122B42800500122B4960122301B0B301A089023CC
:0231330039C082C08332C2AF80FEC2AF301B0B301A0890234DC082C08332C2AF80FE504557
:02315300533A5354455020322E3330004F004500450045004500417272617920496E6465BC
:023173007820546F6F2042696700537472696E67205472756E63617465640046697669641D
:0231930065206279205A65726F004D617468204F766572666C6F77004D61746820556E643B
:0231B3006572666C6F77004D61746820417267756D656E742042616400496E7075742056E1
:0231D300616C7565205466F6F2042696700496E7075742056616C7565204572726F7200491F
:0231F3006E7075742056616C7565204E656761746976650049737461636B2F4235322043F9
:02413006F6E666C6963740049737461636B204572726F7200566563746F722072656C6FEC
:02413300636174696F6E204641494C4544004552524F522000284C696E652000293A200049
:0E245300535441525455502F52455345540096
:00003301CC
```

**Serial snoop processor(IC 16) hex file(Serial5.hex):**

```
:10000000680006061E0A1209080C2A0003042803C9
:100010000606E8051009EA02060A08022900000A95
:10002000540C130A6C0C2B00EB02140A0008010C90
:1000300006000000C0500060C2500000A120C8901C0
:1000400043066504130C890143066505140C8901F8
:1000500043060504150C890143060505160C8901A4
:1000600043062504170C890143062505180C890150
:1000700043064504190C890143064505200C2B0055
:06008000EB02400A000A39
:0203FE00170ADC
:000000001FF
```